



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com

Company: \_\_\_\_\_ Ship To: \_\_\_\_\_  
 Account Number: \_\_\_\_\_ Date: \_\_\_\_\_ Address: \_\_\_\_\_  
 Quote  PO Number: \_\_\_\_\_ City: \_\_\_\_\_ State: \_\_\_\_\_  
 RTS: \_\_\_\_\_ Zip: \_\_\_\_\_ Phone: \_\_\_\_\_  
 Marked For: \_\_\_\_\_ Email/Fax: \_\_\_\_\_

Effective 01/07/2019

## MyWay

**HCPC: E8000**

### Packages

*Includes Frame, Harness, User Handles, Guide Handles, Headrest, Saddle and Ankle Guides. Packages for frame sizes 1 and 2 include Size 1 Saddle; Packages for frame size 3 include Size 2 Saddle.*

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 MyWay Package with Size 1 Harness - Purple	PMWS1P	\$4967
<input type="checkbox"/>	Size 1 MyWay Package with Size 2 Harness - Green	PMWS1G	\$4971
<input type="checkbox"/>	Size 1 MyWay Package with Size 3 Harness - Blue	PMWS1B	\$4999
<input type="checkbox"/>	Size 2 MyWay Package with Size 3 Harness - Blue	PMWS2B	\$5538
<input type="checkbox"/>	Size 2 MyWay Package with Size 4 Harness - Red	PMWS2R	\$5571
<input type="checkbox"/>	Size 2 MyWay Package with Size 5 Harness - Orange	PMWS2O	\$5606
<input type="checkbox"/>	Size 3 MyWay Package with Size 4 Harness - Red	PMWS3R	\$5747
<input type="checkbox"/>	Size 3 MyWay Package with Size 5 Harness - Orange	PMWS3O	\$5782
<input type="checkbox"/>	Size 3 MyWay Package with Size 6 Harness - Gray	PMWS3G	\$5798

### Frame

*Selection required. Includes multi-functional casters.*

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 MyWay Frame	173-1600	\$3488
<input type="checkbox"/>	Size 2 MyWay Frame	173-2600	\$4023
<input type="checkbox"/>	Size 3 MyWay Frame	173-3600	\$4189

### Harness

*Selection required*

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 Harness - Purple	173-4100-04	\$316
<input type="checkbox"/>	Size 2 Harness - Green	173-4200-07	\$320
<input type="checkbox"/>	Size 3 Harness - Blue	173-4300-02	\$348
<input type="checkbox"/>	Size 4 Harness - Red	173-4400-09	\$381
<input type="checkbox"/>	Size 5 Harness - Orange	173-4500-06	\$416
<input type="checkbox"/>	Size 6 Harness - Gray	173-4600-08	\$432

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

## LECKEY<sup>o</sup>

### Handles

*Must match Frame size*

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 User Handles (Pair)	173-1667	\$162
<input type="checkbox"/>	Size 2 User Handles (Pair)	173-2667	\$166
<input type="checkbox"/>	Size 3 User Handles (Pair)	173-3667	\$173
<input type="checkbox"/>	Size 1 Guide Handles	173-1691	\$276
<input type="checkbox"/>	Size 2 Guide Handles	173-2691	\$276
<input type="checkbox"/>	Size 3 Guide Handles	173-3691	\$276

### Positioning Accessories

*Ankle Guides must match Frame size*

	Description	Item Code	Retail
<input type="checkbox"/>	Headrest (One Size)	173-1680	\$197
<input type="checkbox"/>	Size 1 Saddle	173-1671	\$308
<input type="checkbox"/>	Size 2 Saddle	173-2671	\$311
<input type="checkbox"/>	Size 1 Ankle Guides (Pair)	173-1700	\$220
<input type="checkbox"/>	Size 2 Ankle Guides (Pair)	173-2700	\$220
<input type="checkbox"/>	Size 3 Ankle Guides (Pair)	173-3700	\$220

### Other Accessories

*Connection belt available only with Size 2 or Size 3 Frames. Selection of Guide Handles required with Guide Connection Belt.*

	Description	Item Code	Retail
<input type="checkbox"/>	Guide Connection Belt	173-1693	\$100
<input type="checkbox"/>	Hoist Straps (Pair)	173-1800	\$52

# LECKEY<sup>o</sup>

## Technical Specifications

Frame Sizing Specification	Size 1	Size 2	Size 3			
Age (approximate)	1 - 5 years	4 - 11 years	9 - 16 years			
User Shoulder Height (to Floor)	25" - 36"	34" - 47"	42" - 59"			
User Saddle Height (Inside Leg)	8.5" - 17.5"	16" - 24"	21.5" - 32"			
Maximum User Weight	55 lbs	110 lbs	176 lbs			
Frame Weight	26.5 lbs	14.5 lbs	21 lbs			
Upper Support Unit Weight	Size 1 Frame and Upper Support Unit do no separate		23 lbs			
Frame Width	22.5"	26.5"	28.5"			
Frame Length	28.5"	33.5"	43.5"			
Saddle	9"L x 3.25"W	10"L x 3.25"W				
Harness Sizing Specification	Size 1	Size 2	Size 3	Size 4	Size 5	Size 6
Age (approximate)	1 - 4 years	2 - 6 years	5 - 9 years	8 - 12 years	12 - 15 years	16 - 18 years
Waist at Navel Level	18" - 22"	19.5" - 23.5"	21.5" - 25.5"	23.5" - 27.5"	27.5" - 31.5"	31.5" - 35.5"
Color	Purple	Green	Blue	Red	Orange	Gray
Maximum User Weight	110 lbs	110 lbs	110 lbs	176 lbs	176 lbs	176 lbs

## Technical Specifications - Metric

Frame Sizing Specification	Size 1	Size 2	Size 3			
Age (approximate)	1 - 5 years	4 - 11 years	9 - 16 years			
User Shoulder Height (to Floor)	635 mm - 920 mm	870 mm - 1190 mm	1070 mm - 1490 mm			
User Saddle Height (Inside Leg)	220 mm - 440 mm	400 mm - 610 mm	550 mm - 810 mm			
Maximum User Weight	25 kgs	50 kgs	80 kgs			
Frame Weight	12 kgs	6.5 kgs	9.5 kgs			
Upper Support Unit Weight	Size 1 Frame and Upper Support Unit do no separate		10.5 kgs			
Frame Width	570 mm	670 mm	720 mm			
Frame Length	720 mm	850 mm	1100 mm			
Saddle	225 mm L x 80 mm W	255 mm L x 80 mm W				
Harness Sizing Specification	Size 1	Size 2	Size 3	Size 4	Size 5	Size 6
Age (approximate)	1 - 4 years	2 - 6 years	5 - 9 years	8 - 12 years	12 - 15 years	16 - 18 years
Waist at Navel Level	460 mm - 560 mm	500 mm - 600 mm	550 mm - 650 mm	600 mm - 700 mm	700 mm - 800 mm	800 mm - 900 mm
Color	Purple	Green	Blue	Red	Orange	Gray
Maximum User Weight	50 kgs	50 kgs	50 kgs	80 kgs	80 kgs	80 kgs