



Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com

Company: _____

Ship To: _____

Account Number: _____ Date: _____

Address: _____

Quote PO Number: _____

City: _____ State: _____

RTS: _____

Zip: _____ Phone: _____

Marked For: _____

Email/Fax: _____

Effective 08/27/2018

Mygo Seating System

Size 1 Weight Capacity: 110 lbs | Size 2 Weight Capacity: 132 lbs

Seat Shell

Selection required

	Description	Item Code	HCPC	Retail
	Size 1 Mygo Seat Shell	117-610		
	Seat Base with Depth Adjustment	117-610-1	E2294	\$768
	Seat Back - Fixed	117-610-2	E2293	\$768
<input type="checkbox"/>	Seat Back Reclining Mechanism	117-610-3	E1014	\$404
	Seat Base Hip Guides/Support	117-610-4	E0956	\$112
	Seat Base Medial Thigh Guides	117-610-5	E0957	\$139
	Seat Base Leg Adductors	117-610-6	K0108	\$112
	Size 2 Mygo Seat Shell	137-600		
	Seat Base with Depth Adjustment	137-600-1	E2294	\$851
	Seat Back - Fixed	137-600-2	E2293	\$851
<input type="checkbox"/>	Seat Back Reclining Mechanism	137-600-3	E1014	\$404
	Seat Base Hip Guides/Support	137-600-4	E0956	\$112
	Seat Base Medial Thigh Guides	137-600-5	E0957	\$139
	Seat Base Leg Adductors	137-600-6	K0108	\$112

Seat Cover

Selection required. Must match Seat Shell size.

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Size 1 Mygo Seat Cover - Blue	117-1851-08		N/C
<input type="checkbox"/>	Size 1 Mygo Seat Cover - Gray	117-1851-06		N/C
<input type="checkbox"/>	Size 1 Mygo Seat Cover - Orange	117-1851-07		N/C
<input type="checkbox"/>	Size 1 Mygo Seat Cover - Pink	117-1851-09		N/C
<input type="checkbox"/>	Size 2 Mygo Seat Cover - Blue	117-2851-08		N/C
<input type="checkbox"/>	Size 2 Mygo Seat Cover - Gray	117-2851-06		N/C
<input type="checkbox"/>	Size 2 Mygo Seat Cover - Orange	117-2851-07		N/C
<input type="checkbox"/>	Size 2 Mygo Seat Cover - Pink	117-2851-09		N/C

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

LECKEY^o

Base

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Size 2 Hi Low Chassis <i>Foot Pedal Operated. Includes Push Handle.</i>	102-2753	E1399	\$1895
<input type="checkbox"/>	Size 2 Kimba Neo - Black <i>Possible interference with size 2 Mygo Seat</i>	470G7120000350	E1234	\$2205
	Transit System for Kimba	471X00ST007	E1399	\$51

Pelvic Fixation

Selection of either Pelvic Harness or Pelvic Cradle required

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Pelvic Harness & Brackets Small - Black	117-801	E0978	\$137
<input type="checkbox"/>	Pelvic Harness & Brackets Medium - Black	117-802	E0978	\$137
<input type="checkbox"/>	Pelvic Harness Spacer Pads - Black (Pair)	120-789	K0108	\$78
<input type="checkbox"/>	Pelvic Cradle Extra Small - Black	133-727	E0978	\$340
<input type="checkbox"/>	Pelvic Cradle Small - Black	133-728	K0108	\$383
<input type="checkbox"/>	Pelvic Cradle Medium - Black	133-729	K0108	\$425
<input type="checkbox"/>	Hip Laterals (Pair) <i>Pelvic Cradle required</i>	137-632	E0956	\$284
<input type="checkbox"/>	Obliquity Kit	117-791	K0108	\$139

Dynamic Backrest Kit

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Dynamic Backrest 50N Cylinder	137-660	K0108	\$391
<input type="checkbox"/>	Dynamic Backrest 100N Cylinder	137-661	K0108	\$391

Armrest

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Size 1 Mygo Adjustable Armrests (Pair)	117-1777	E0973	\$280
<input type="checkbox"/>	Size 2 Mygo Adjustable Armrests (Pair)	117-2777	E0973	\$308

LECKEY^o

Lower Leg Positioning

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Size 1 Mygo Angled Footplates (Pair)	117-1776	K0040	\$280
<input type="checkbox"/>	Size 2 Mygo Angled Footplates (Pair)	117-2776	K0040	\$308
<input type="checkbox"/>	Small Sandals with Straps (Pair)	152-1600	E0951, E0952	\$225
<input type="checkbox"/>	Medium Sandals with Straps (Pair)	152-2600	E0951, E0952	\$245
<input type="checkbox"/>	Large Sandals with Straps (Pair)	152-3600	E0951, E0952	\$265
<input type="checkbox"/>	Sandal Raiser 1" (Each)	Quantity: _____ (1-6)	117-787	\$76
<input type="checkbox"/>	Knee Pads (Pair)		117-765	\$139

Harness

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Small Trunk Harness	120-799	E0960	\$249
<input type="checkbox"/>	Medium Trunk Harness	120-757	E0960	\$249
<input type="checkbox"/>	Large Trunk Harness	117-761	E0960	\$249

Lateral Support

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Flip-Away Laterals with Cover (Pair)	117-837	E0956, E1028	\$488
<input type="checkbox"/>	Flip-Away Laterals Extended Height with Cover (Pair)	137-634	E0956, E1028	\$599
<input type="checkbox"/>	Rigid Laterals with Cover (Pair)	117-803	E0956	\$279
<input type="checkbox"/>	Rigid Laterals Extended Height with Cover (Pair)	137-633	E0956	\$390
<input type="checkbox"/>	Flip-Away Complex Lateral Horizontal with Cover Left	137-637	E0956, E1028	\$312
<input type="checkbox"/>	Flip-Away Complex Lateral Horizontal with Cover Right	137-638	E0956, E1028	\$312
<input type="checkbox"/>	Flip-Away Complex Lateral Vertical with Cover Left	137-635	E0956, E1028	\$312
<input type="checkbox"/>	Flip-Away Complex Lateral Vertical with Cover Right	137-636	E0956, E1028	\$312
<input type="checkbox"/>	Chest Support - Black <i>Requires laterals. Not available with Complex Laterals.</i>	137-808-05	E0960	\$157

Headrest

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Contoured Headrest Support Hardware	120-859	E1028	\$317
	Contoured Headrest Cushion - Black	137-800-05	E0955	\$145
<input type="checkbox"/>	Flat Headrest Support Hardware	117-697	E1028	\$194
	Flat Headrest Cushion - Black	117-850-05	E0955	\$98
<input type="checkbox"/>	Flat Headrest Laterals with Covers - Black (Pair) <i>Requires Flat Headrest. Not available with Contoured Headrest.</i>	120-893-05	E0956	\$278
<input type="checkbox"/>	Removable Head Support Bracket for 1/2" Box Section	117-745	E1028	\$62

LECKEY^o

Tray

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Tray - Gray	117-724	E0950	\$335
<input type="checkbox"/>	Padded Tray Insert - Gray	117-270	K0108	\$291
<input type="checkbox"/>	Grab Rail/Toy Bar for Tray	117-769	K0108	\$251

Kimba Neo Accessories

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Kimba Neo Sun/Rain Canopy Complete Gray/Black <i>Kimba Neo Canopy Hardware required</i>	471S71SK301	K0108	\$403
<input type="checkbox"/>	Kimba Neo Size 2 Canopy Hardware	415S88SK024		\$84
<input type="checkbox"/>	Kimba Neo Basic Canopy Black <i>Kimba Neo Canopy Hardware required</i>	471S71SK505	K0108	\$175
<input type="checkbox"/>	Kimba Neo Size 2 Canopy Hardware	415S88SK024		\$84
<input type="checkbox"/>	Kimba Neo Vent Mount for Pulmonetics LTV	471C49NTKNLTVM	E1029/E1028	\$545
<input type="checkbox"/>	Kimba Neo Vent Mount for Respironics Trilogy	471C49NTKNWMTRI	E1029/E1028	\$545
<input type="checkbox"/>	Kimba Neo Size 2 Horizontal Vent Tray	471S71SK160	E1029	\$464
<input type="checkbox"/>	Kimba Neo IV Pole	471C49NTKNIVP48	E0776	\$423
<input type="checkbox"/>	Kimba Neo O2 Tank Holder Side Mount	471C49NTKNO2H	E2208	\$261

Accessories

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	1" Back Spacer Pad and Rigid Lateral Support - Blue	117-830-07	E0956, E1297	\$540
<input type="checkbox"/>	1" Back Spacer Pad and Rigid Lateral Support - Gray	117-830-09	E0956, E1297	\$540
<input type="checkbox"/>	1" Back Spacer Pad and Rigid Lateral Support - Orange	117-830-08	E0956, E1297	\$540
<input type="checkbox"/>	1" Back Spacer Pad and Rigid Lateral Support - Pink	117-830-06	E0956, E1297	\$540
<input type="checkbox"/>	Shoulder Support Hardware for Mygo Size 1 <i>Shoulder Support Cushion required. Shoulder Support standard on Size 2.</i>	137-621	K0108	\$237
<input type="checkbox"/>	Shoulder Support Cushion - Blue	137-801-08	K0108	\$135
<input type="checkbox"/>	Shoulder Support Cushion - Gray	137-801-06	K0108	\$135
<input type="checkbox"/>	Shoulder Support Cushion - Orange	137-801-07	K0108	\$135
<input type="checkbox"/>	Shoulder Support Cushion - Pink	137-801-09	K0108	\$135

Assembly

Product Assembly Fee not subject to discounts

<input type="checkbox"/>	Product Assembly Fee			\$120
--------------------------	----------------------	--	--	-------

LECKEY^o

Universal Specifications

Specification	Size 1	Size 2
Age (approximate)	3 - 10 years	8 - 14 years
Backrest Angle: Prone	10°	10°
Backrest Angle: Recline	25°	25°
Footplate: Abduction	8°	8°
Footplate: Adduction	12°	20°
Footplate: Plantar/Dorsiflexion	10° / 10°	10° / 10°
Size 2 Hi Lo Chassis (Excluding Push Bar)		
Tilt in Space: Prone	10°	10°
Tilt in Space: Recline	25°	25°

Technical Specifications

Specification	Size 1	Size 2
User Height	41.3" - 59"	50" - 66"
User Weight	39.6 - 110 lbs	39.6 - 132 lbs
Seat Depth	10.6" - 16.5"	13.8" - 18.5"
Seat Width	8" - 12.8"	8.7" - 13.6"
Backrest Height	14.2" - 18.5"	18.1" - 22.4"
Seat to Sandal	8.5" - 13.8"	12.4" - 18.5"
Chest Width (Distance between laterals)	6.7" - 10.6"	6.7" - 10.6"
Armrest Height	6.3" - 8.3"	8.3" - 10.2"
Knee Width	3.5" - 4.3"	4.7" - 5.5"
Tray Size	21.6" x 18.9"	21.6" x 18.9"
Seat Unit Weight	22 lbs	31.9 lbs
Size 2 Hi Lo Chassis (Excluding Push Bar)		
Top of Seat to Floor	13.4" - 25.8"	13.4" - 25.8"
Base Weight	26.5 lbs	26.5 lbs

Technical Specifications - Metric

Specification	Size 1	Size 2
User Height	1050 mm - 1500 mm	1270 mm - 1680 mm
User Weight	18 - 50 kg	18 - 60 kg
Seat Depth	270 mm - 420 mm	350 mm - 470 mm
Seat Width	200 mm - 325 mm	220 mm - 345 mm
Backrest Height	360 mm - 470 mm	460 mm - 570 mm
Seat to Sandal	215 mm - 350 mm	315 mm - 470 mm
Chest Width (Distance between laterals)	170 mm - 270 mm	170 mm - 270 mm
Armrest Height	160 mm - 210 mm	210 mm - 260 mm
Knee Width	90 mm - 110 mm	120 mm - 140 mm
Tray Size	550 mm x 480 mm	550 mm x 480 mm
Seat Unit Weight	10 kg	14.5 kg
Size 2 Hi Lo Chassis (Excluding Push Bar)		
Top of Seat to Floor	340 mm - 655 mm	340 mm - 655 mm
Base Weight	12 kg	12 kg