



Pronestander

Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com

Company: _____ Ship To: _____
 Account Number: _____ Date: _____ Address: _____
 Quote PO Number: _____ City: _____ State: _____
 RTS: _____ Zip: _____ Phone: _____
 Marked For: _____ Email/Fax: _____

Effective 08/27/2018

Pronestander

HCPC: E0638

Standard Model

Selection Required. Includes Pelvic Band, Flexible Laterals, PU Knee Supports, Footplate, Sandals, Tray, Tray Cover and Activity Bowl.

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 Pronestander - Meadow	LPS/1-MW-EXP	\$2675
<input type="checkbox"/>	Size 2 Pronestander - Mondeo Blue	LPS/2-MB-EXP	\$2925
<input type="checkbox"/>	Size 3 Pronestander - Mondeo Blue	LPS/3-MB-EXP	\$3425

Headrest

	Description	Item Code	Retail
<input type="checkbox"/>	Headrest Interface Frame <i>Required with Contoured Headrest</i>	LIF/1	\$470
<input type="checkbox"/>	Size 1 Contoured Headrest - Vinyl	LIF/2/A	\$319
<input type="checkbox"/>	Size 1 Contoured Headrest - Meadow	LIF/2/A-MW	\$319
<input type="checkbox"/>	Size 2 Contoured Headrest - Vinyl	LIF/2/B	\$319
<input type="checkbox"/>	Size 2 Contoured Headrest - Mondeo Blue	LIF/2/B-MB	\$319
<input type="checkbox"/>	Size 1 Hyper Extension Pad - Mondeo Blue	LIF/3/A-MB	\$355
<input type="checkbox"/>	Size 2 Hyper Extension Pad - Mondeo Blue	LIF/3/B-MB	\$405

Additional Pelvic Band

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 Pelvic Band - Meadow	LPS/1/B-MW	\$160
<input type="checkbox"/>	Size 2 Pelvic Band - Mondeo Blue	LPS/2/B-MB	\$175
<input type="checkbox"/>	Size 3 Pelvic Band - Mondeo Blue	LPS/3/B-MB	\$190

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

LECKEY^o

Additional Knee Supports

Description	Item Code	Retail (Each)
<input type="checkbox"/> Size 1 Upholstered Knee Cups - Meadow (Pair)	LPS/1/K-MW	\$109
<input type="checkbox"/> Size 2 Upholstered Knee Cups - Mondeo Blue (Pair)	LPS/2/K-MB	\$109
<input type="checkbox"/> Size 3 Upholstered Knee Cups - Mondeo Blue (Pair)	LPS/3/K-MB	\$109
<input type="checkbox"/> Size 4 Upholstered Knee Cups - Mondeo Blue (Pair)	LPS/4/K-MB	\$109
<input type="checkbox"/> Size 1 Knee Spacer (Each)	<input type="checkbox"/> Quantity 1 <input type="checkbox"/> Quantity 2	LPS/1/F \$55
<input type="checkbox"/> Size 2 Knee Spacer (Each)	<input type="checkbox"/> Quantity 1 <input type="checkbox"/> Quantity 2	LPS/2/F \$55
<input type="checkbox"/> Size 3 Knee Spacer (Each)	<input type="checkbox"/> Quantity 1 <input type="checkbox"/> Quantity 2	LPS/3/F \$55

Additional Foot Positioning

Description	Item Code	Retail
<input type="checkbox"/> Small Sandals with Straps (Pair)	152-1600	\$225
<input type="checkbox"/> Medium Sandals with Straps (Pair)	152-2600	\$245
<input type="checkbox"/> Large Sandals with Straps (Pair)	152-3600	\$265
<input type="checkbox"/> Size 1 Pronestander Sandal Raiser (Each)	Quantity: _____ (1-6)	LPS/1/J \$62
<input type="checkbox"/> Size 2 Pronestander Sandal Raiser (Each)	Quantity: _____ (1-6)	LPS/2/J \$62
<input type="checkbox"/> Size 3 Pronestander Sandal Raiser (Each)	Quantity: _____ (1-6)	LPS/3/J \$62

Accessories

Description	Item Code	Retail
<input type="checkbox"/> Size 1 Grab/Suspension Rail	LPS/1/C	\$179
<input type="checkbox"/> Size 2 Grab/Suspension Rail	LPS/2/C	\$179
<input type="checkbox"/> Size 3 Grab/Suspension Rail	LPS/3/C	\$179
<input type="checkbox"/> Size 1 Pronestander Grab Post	LPS/1/D	\$179
<input type="checkbox"/> Size 2 Pronestander Grab Post	LPS/2/D	\$179
<input type="checkbox"/> Size 3 Pronestander Grab Post	LPS/3/D	\$179

LECKEY[®]

Technical Specifications

Specification	Size 1	Size 2	Size 3
Age (approximate)	1 - 5 years	4 - 10 years	9 - 18 years
User Height (Min / Max)	31.5" / 45.3"	43.3" / 57"	55" / 73"
Maximum User Weight	132 lbs	176 lbs	198 lbs
Chest Plate Height	18.5" - 27"	27" - 38"	38" - 53"
Base	32" x 24"	42" x 26"	50" x 28"
Activity Tray	16" x 21"	18" x 23"	20" x 25"

Technical Specifications - Metric

Specification	Size 1	Size 2	Size 3
Age (approximate)	1 - 5 years	4 - 10 years	9 - 18 years
User Height (Min / Max)	800 mm / 1150 mm	1100 mm / 1450 mm	1400 mm / 1850 mm
Maximum User Weight	60 kg	80 kg	90 kg
Chest Plate Height	470 mm - 680 mm	680 mm - 980 mm	970 mm - 1140 mm
Base	820 mm x 600 mm	1070 mm x 650 mm	1270 mm x 700 mm
Activity Tray	400 mm x 525 mm	450 mm x 575 mm	500 mm x 600 mm