



Phone (800) 981-1540 | Fax (715) 254-0996  
sales@kimobility.com | www.kimobility.com

# ROGUE DT

Company:

Account Number:  Date:

☐ Quote: ☐ PO Number:

ATP:

Marked For:

Associated Quote/Order#:

☐ Rogue DT without Transit (K0005, E1235) \$2930

☐ Rogue DT with Transit (K0005, E1235) \$3277

300 lb Weight Capacity

Ship To:

Address:

City:  State:

Zip:  Phone:

Email/Fax:

**Submit Form**

Effective: 8/17/2026

## Power Assist

☐ Twist Drive System (E0986) \$7950

*Includes Multi-Function Toggle Controller, Removable Li-ion Battery - 36V, 4Ah, Twist Carry Bag, and Rear Mounting Hardware. See page 6 for additional, required selections.  
Weight Capacity: Front - 265lb (120 kg), Rear - 309lb (140 kg).*

## Frame

### Seat Width

*See reference page (A). 12" and 13" Seat Widths not available with Klaxon Twist.*

- |                                      |       |
|--------------------------------------|-------|
| <input type="checkbox"/> 12"         | N/C   |
| <input type="checkbox"/> 13"         | N/C   |
| <input type="checkbox"/> 14"         | N/C   |
| <input type="checkbox"/> 15"         | N/C   |
| <input type="checkbox"/> 16"         | N/C   |
| <input type="checkbox"/> 17"         | N/C   |
| <input type="checkbox"/> 18"         | N/C   |
| <input type="checkbox"/> 19"         | N/C   |
| <input type="checkbox"/> 20" (E2201) | \$575 |

### Seat Depth

*See reference page (B)*

- |                                      |       |
|--------------------------------------|-------|
| <input type="checkbox"/> 12"         | N/C   |
| <input type="checkbox"/> 13"         | N/C   |
| <input type="checkbox"/> 14"         | N/C   |
| <input type="checkbox"/> 15"         | N/C   |
| <input type="checkbox"/> 16"         | N/C   |
| <input type="checkbox"/> 17"         | N/C   |
| <input type="checkbox"/> 18"         | N/C   |
| <input type="checkbox"/> 19"         | N/C   |
| <input type="checkbox"/> 20" (E2203) | \$575 |

### Transit Option

*Not retrofittable. Transit not available with Non-Folding Backrest.*

- |                                                 |       |
|-------------------------------------------------|-------|
| <input type="checkbox"/> Transit System (E1022) | \$347 |
| <input type="checkbox"/> Omit                   | N/C   |

### Seat Option

*Solid Seat Pan not available with Ergonomic Seat Depth*

- |                                                          |       |
|----------------------------------------------------------|-------|
| <input type="checkbox"/> Standard Seat Sling             | N/C   |
| <input type="checkbox"/> Aluminum Solid Seat Pan (E2231) | \$277 |

### Adjustable Depth

*No Rear Adjustment selection required with Non-Folding Backrest*

- |                                             |       |
|---------------------------------------------|-------|
| <input type="checkbox"/> No Rear Adjustment | N/C   |
| <input type="checkbox"/> 1/2"               | N/C   |
| <input type="checkbox"/> 1"                 | \$89  |
| <input type="checkbox"/> 2"                 | \$172 |

### Sling Position

*See reference page (C)*

- |                                           |     |
|-------------------------------------------|-----|
| <input type="checkbox"/> 0" (Performance) | N/C |
| <input type="checkbox"/> 1"               | N/C |
| <input type="checkbox"/> 2" (Classic)     | N/C |

### Front Frame Bend

*See reference page (D). 6" Casters not available with 85° & 90°.*

- |                              |     |
|------------------------------|-----|
| <input type="checkbox"/> 90° | N/C |
| <input type="checkbox"/> 85° | N/C |
| <input type="checkbox"/> 80° | N/C |
| <input type="checkbox"/> 75° | N/C |

### Ergonomic Seat Depth

*See reference page (E). Recommended to use with air cushions. Ergonomic Seat is N/C.*

- |                               |     |
|-------------------------------|-----|
| <input type="checkbox"/> Omit | N/C |
|-------------------------------|-----|

|            |    | Use Seat Depth Selection |                          |                          |                          |                          |                          |                          |                          |                          |
|------------|----|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|
|            |    | 12"                      | 13"                      | 14"                      | 15"                      | 16"                      | 17"                      | 18"                      | 19"                      | 20"                      |
| Ergo Depth | 5" | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
|            | 6" | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
|            | 7" |                          | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
|            | 8" |                          |                          | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
|            | 9" |                          |                          |                          | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |

HCPCS Codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS Codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN1408.1 Pricing and specifications subject to change without notice.



Phone (800) 981-1540 | Fax (715) 254-0996  
sales@kimobility.com | www.kimobility.com

# ROGUE DT

## Backrest

### Backrest Type

**Non-Folding Backrest required for Bolt-On and Fixed Side Guard Mount. Non-Folding Backrest not available with Transit or Swing Away Armrests. Non-Folding Backrest in combination with Fixed Height Back Post, maximum 16" Back Height. Non-Folding Backrest is 1.3 lbs lighter.**

- ☐ Folding/Locking Backrest \_\_\_\_\_ N/C  
☐ Non-Folding Backrest (Welded) \_\_\_\_\_ \$110

### Back Post / Height Range / Back Height

**Fixed Height Back Post not available with Back Upholstery with Folding/Locking Backrest. Fold Down Push Handle not available with Back Heights 11" and under.**

| Back Post         | Price | Height Range                    | Back Height                  |                              |                              |                              |                              |
|-------------------|-------|---------------------------------|------------------------------|------------------------------|------------------------------|------------------------------|------------------------------|
| Height Adjustable | N/C   | <input type="checkbox"/> Short  | <input type="checkbox"/> 11" | <input type="checkbox"/> 12" | <input type="checkbox"/> 13" | <input type="checkbox"/> 14" |                              |
|                   |       | <input type="checkbox"/> Medium | <input type="checkbox"/> 13" | <input type="checkbox"/> 14" | <input type="checkbox"/> 15" | <input type="checkbox"/> 16" | <input type="checkbox"/> 17" |
|                   |       | <input type="checkbox"/> Tall   | <input type="checkbox"/> 17" | <input type="checkbox"/> 18" | <input type="checkbox"/> 19" | <input type="checkbox"/> 20" |                              |
| Fixed Height      | N/C   | <input type="checkbox"/> Fixed  | <input type="checkbox"/> 9"  | <input type="checkbox"/> 10" | <input type="checkbox"/> 11" | <input type="checkbox"/> 12" | <input type="checkbox"/> 13" |
|                   |       |                                 | <input type="checkbox"/> 14" | <input type="checkbox"/> 15" | <input type="checkbox"/> 16" | <input type="checkbox"/> 17" | <input type="checkbox"/> 18" |

### Seat Back Angle

**-6°, +/-9°, -12° not available with Non-Folding Backrest**

- ☐ Standard \_\_\_\_\_ N/C  
☐ +3° \_\_\_\_\_ N/C  
☐ +6° \_\_\_\_\_ N/C  
☐ +9° \_\_\_\_\_ N/C  
☐ -3° \_\_\_\_\_ N/C  
☐ -6° \_\_\_\_\_ N/C  
☐ -9° \_\_\_\_\_ N/C  
☐ -12° \_\_\_\_\_ N/C

### Back Upholstery

**Not available with Fixed Height Back Post with Folding/Locking Backrest or 9" - 10" Back Heights. If both Upholstery and Axiom Back are selected, the Upholstery will be installed.**

- ☐ Padded Velcro Adjustable with Standard Cover \_\_\_\_\_ N/C  
☐ Padded Velcro Adjustable with Air Mesh Breathable Cover \_\_\_\_\_ \$81  
☐ Tension Adjustable with Standard Cover (E2611) \_\_\_\_\_ \$374  
☐ Tension Adjustable with Air Mesh Breathable Cover (E2611) \_\_\_\_\_ \$416  
☐ Omit \_\_\_\_\_ (\$10)

### Push Handles

**Low Profile Push Handle is 2" shorter than Standard Push Handle. Bolt-On Push Handle will be 1" - 2" above Back Height.**

- ☐ Omit Push Handle \_\_\_\_\_ N/C  
☐ Standard Push Handle \_\_\_\_\_ N/C  
☐ Low Profile Push Handle \_\_\_\_\_ \$157  
☐ Bolt-On Push Handle \_\_\_\_\_ \$183  
☐ Fold Down Push Handle \_\_\_\_\_ \$115

### Stroller Handle

**Stroller Handles are removable and height adjustable**

- ☐ Single Stroller Handle \_\_\_\_\_ \$172  
☐ Double Stroller Handle \_\_\_\_\_ \$345

### Rigidizer Bar

**Adjustable Height not available with Non-Folding Backrest. Non-Adjustable Height not available with Fixed Height Back Post with Folding/Locking Backrest.**

- ☐ Adjustable Height Rigidizer Bar \_\_\_\_\_ \$79  
☐ Non-Adjustable Height Rigidizer Bar \_\_\_\_\_ N/C

### Axiom Back - Easy Fit

**Optimized by Ki to fit the Seat Width and Back Height selection. See Axiom Order Form if you require different sizes and selections. Pricing reflects only the Easy Fit option, standard pricing still applies for Axiom Backs.**

### Back Type

- ☐ Axiom Posterior Back \_\_\_\_\_ N/C

### Back Hardware

- ☐ DEX Hardware \_\_\_\_\_ N/C  
☐ LWX Hardware \_\_\_\_\_ N/C

### Military Service Patches

- ☐ Army ☐ Navy ☐ Marines ☐ Air Force  
☐ Coast Guard ☐ US Purple Heart \_\_\_\_\_ \$94  
☐ Combo Army/Purple Heart ☐ Combo Navy/Purple Heart  
☐ Combo Marines/Purple Heart ☐ Combo Air Force/Purple Heart  
☐ Combo Coast Guard/Purple Heart \_\_\_\_\_ \$152

## Front Seat Height and Casters

### Caster Size / Footrest Taper / Front Seat Height

See reference page (G, I). Select Taper, Caster, then Front Seat Height. 2" Taper or V-Taper not available with 75° Front Frame Bend.

| Footrest Taper                        | Price | Caster Size | Front Seat Height              |                                |                                |                                |                                |                                |                                |                                |                                |                                |                                |                                |                                |                                |                                |  |
|---------------------------------------|-------|-------------|--------------------------------|--------------------------------|--------------------------------|--------------------------------|--------------------------------|--------------------------------|--------------------------------|--------------------------------|--------------------------------|--------------------------------|--------------------------------|--------------------------------|--------------------------------|--------------------------------|--------------------------------|--|
| <input type="checkbox"/> Straight     | N/C   | 3"          | <input type="checkbox"/> 13.5" | <input type="checkbox"/> 14"   | <input type="checkbox"/> 14.5" | <input type="checkbox"/> 15"   | <input type="checkbox"/> 15.5" | <input type="checkbox"/> 16"   | <input type="checkbox"/> 16.5" | <input type="checkbox"/> 17"   | <input type="checkbox"/> 17.5" | <input type="checkbox"/> 18"   | <input type="checkbox"/> 18.5" | <input type="checkbox"/> 19"   | <input type="checkbox"/> 19.5" | <input type="checkbox"/> 20"   | <input type="checkbox"/> 20.5" |  |
|                                       |       | 4"          | <input type="checkbox"/> 14"   | <input type="checkbox"/> 14.5" | <input type="checkbox"/> 15"   | <input type="checkbox"/> 15.5" | <input type="checkbox"/> 16"   | <input type="checkbox"/> 16.5" | <input type="checkbox"/> 17"   | <input type="checkbox"/> 17.5" | <input type="checkbox"/> 18"   | <input type="checkbox"/> 18.5" | <input type="checkbox"/> 19"   | <input type="checkbox"/> 19.5" | <input type="checkbox"/> 20"   | <input type="checkbox"/> 20.5" |                                |  |
|                                       |       | 5"          | <input type="checkbox"/> 14.5" | <input type="checkbox"/> 15"   | <input type="checkbox"/> 15.5" | <input type="checkbox"/> 16"   | <input type="checkbox"/> 16.5" | <input type="checkbox"/> 17"   | <input type="checkbox"/> 17.5" | <input type="checkbox"/> 18"   | <input type="checkbox"/> 18.5" | <input type="checkbox"/> 19"   | <input type="checkbox"/> 19.5" | <input type="checkbox"/> 20"   | <input type="checkbox"/> 20.5" |                                |                                |  |
|                                       |       | 6"          | <input type="checkbox"/> 15"   | <input type="checkbox"/> 15.5" | <input type="checkbox"/> 16"   | <input type="checkbox"/> 16.5" | <input type="checkbox"/> 17"   | <input type="checkbox"/> 17.5" | <input type="checkbox"/> 18"   | <input type="checkbox"/> 18.5" | <input type="checkbox"/> 19"   | <input type="checkbox"/> 19.5" | <input type="checkbox"/> 20"   | <input type="checkbox"/> 20.5" |                                |                                |                                |  |
| <input type="checkbox"/> 1"           | N/C   | 3"          | <input type="checkbox"/> 16.5" | <input type="checkbox"/> 17"   | <input type="checkbox"/> 17.5" | <input type="checkbox"/> 18"   | <input type="checkbox"/> 18.5" | <input type="checkbox"/> 19"   | <input type="checkbox"/> 19.5" | <input type="checkbox"/> 20"   | <input type="checkbox"/> 20.5" |                                |                                |                                |                                |                                |                                |  |
|                                       |       | 4"          | <input type="checkbox"/> 17"   | <input type="checkbox"/> 17.5" | <input type="checkbox"/> 18"   | <input type="checkbox"/> 18.5" | <input type="checkbox"/> 19"   | <input type="checkbox"/> 19.5" | <input type="checkbox"/> 20"   | <input type="checkbox"/> 20.5" |                                |                                |                                |                                |                                |                                |                                |  |
|                                       |       | 5"          | <input type="checkbox"/> 17.5" | <input type="checkbox"/> 18"   | <input type="checkbox"/> 18.5" | <input type="checkbox"/> 19"   | <input type="checkbox"/> 19.5" | <input type="checkbox"/> 20"   | <input type="checkbox"/> 20.5" |                                |                                |                                |                                |                                |                                |                                |                                |  |
|                                       |       | 6"          | <input type="checkbox"/> 19"   | <input type="checkbox"/> 19.5" | <input type="checkbox"/> 20"   | <input type="checkbox"/> 20.5" |                                |                                |                                |                                |                                |                                |                                |                                |                                |                                |                                |  |
| <input type="checkbox"/> 2"           | N/C   | 3"          | <input type="checkbox"/> 16.5" | <input type="checkbox"/> 17"   | <input type="checkbox"/> 17.5" | <input type="checkbox"/> 18"   | <input type="checkbox"/> 18.5" | <input type="checkbox"/> 19"   | <input type="checkbox"/> 19.5" | <input type="checkbox"/> 20"   | <input type="checkbox"/> 20.5" |                                |                                |                                |                                |                                |                                |  |
|                                       |       | 4"          | <input type="checkbox"/> 17"   | <input type="checkbox"/> 17.5" | <input type="checkbox"/> 18"   | <input type="checkbox"/> 18.5" | <input type="checkbox"/> 19"   | <input type="checkbox"/> 19.5" | <input type="checkbox"/> 20"   | <input type="checkbox"/> 20.5" |                                |                                |                                |                                |                                |                                |                                |  |
|                                       |       | 5"          | <input type="checkbox"/> 17.5" | <input type="checkbox"/> 18"   | <input type="checkbox"/> 18.5" | <input type="checkbox"/> 19"   | <input type="checkbox"/> 19.5" | <input type="checkbox"/> 20"   | <input type="checkbox"/> 20.5" |                                |                                |                                |                                |                                |                                |                                |                                |  |
|                                       |       | 6"          | <input type="checkbox"/> 19"   | <input type="checkbox"/> 19.5" | <input type="checkbox"/> 20"   | <input type="checkbox"/> 20.5" |                                |                                |                                |                                |                                |                                |                                |                                |                                |                                |                                |  |
| <input type="checkbox"/> 3.5° V-Front | N/C   | 3"          | <input type="checkbox"/> 13.5" | <input type="checkbox"/> 14"   | <input type="checkbox"/> 14.5" | <input type="checkbox"/> 15"   | <input type="checkbox"/> 15.5" | <input type="checkbox"/> 16"   | <input type="checkbox"/> 16.5" | <input type="checkbox"/> 17"   | <input type="checkbox"/> 17.5" | <input type="checkbox"/> 18"   | <input type="checkbox"/> 18.5" | <input type="checkbox"/> 19"   | <input type="checkbox"/> 19.5" | <input type="checkbox"/> 20"   | <input type="checkbox"/> 20.5" |  |
|                                       |       | 4"          | <input type="checkbox"/> 14"   | <input type="checkbox"/> 14.5" | <input type="checkbox"/> 15"   | <input type="checkbox"/> 15.5" | <input type="checkbox"/> 16"   | <input type="checkbox"/> 16.5" | <input type="checkbox"/> 17"   | <input type="checkbox"/> 17.5" | <input type="checkbox"/> 18"   | <input type="checkbox"/> 18.5" | <input type="checkbox"/> 19"   | <input type="checkbox"/> 19.5" | <input type="checkbox"/> 20"   | <input type="checkbox"/> 20.5" |                                |  |
|                                       |       | 5"          | <input type="checkbox"/> 14.5" | <input type="checkbox"/> 15"   | <input type="checkbox"/> 15.5" | <input type="checkbox"/> 16"   | <input type="checkbox"/> 16.5" | <input type="checkbox"/> 17"   | <input type="checkbox"/> 17.5" | <input type="checkbox"/> 18"   | <input type="checkbox"/> 18.5" | <input type="checkbox"/> 19"   | <input type="checkbox"/> 19.5" | <input type="checkbox"/> 20"   | <input type="checkbox"/> 20.5" |                                |                                |  |
|                                       |       | 6"          | <input type="checkbox"/> 15"   | <input type="checkbox"/> 15.5" | <input type="checkbox"/> 16"   | <input type="checkbox"/> 16.5" | <input type="checkbox"/> 17"   | <input type="checkbox"/> 17.5" | <input type="checkbox"/> 18"   | <input type="checkbox"/> 18.5" | <input type="checkbox"/> 19"   | <input type="checkbox"/> 19.5" | <input type="checkbox"/> 20"   | <input type="checkbox"/> 20.5" |                                |                                |                                |  |
| <input type="checkbox"/> 7° V-Front   | N/C   | 3"          | <input type="checkbox"/> 13.5" | <input type="checkbox"/> 14"   | <input type="checkbox"/> 14.5" | <input type="checkbox"/> 15"   | <input type="checkbox"/> 15.5" | <input type="checkbox"/> 16"   | <input type="checkbox"/> 16.5" | <input type="checkbox"/> 17"   | <input type="checkbox"/> 17.5" | <input type="checkbox"/> 18"   | <input type="checkbox"/> 18.5" | <input type="checkbox"/> 19"   | <input type="checkbox"/> 19.5" | <input type="checkbox"/> 20"   | <input type="checkbox"/> 20.5" |  |
|                                       |       | 4"          | <input type="checkbox"/> 14"   | <input type="checkbox"/> 14.5" | <input type="checkbox"/> 15"   | <input type="checkbox"/> 15.5" | <input type="checkbox"/> 16"   | <input type="checkbox"/> 16.5" | <input type="checkbox"/> 17"   | <input type="checkbox"/> 17.5" | <input type="checkbox"/> 18"   | <input type="checkbox"/> 18.5" | <input type="checkbox"/> 19"   | <input type="checkbox"/> 19.5" | <input type="checkbox"/> 20"   | <input type="checkbox"/> 20.5" |                                |  |
|                                       |       | 5"          | <input type="checkbox"/> 14.5" | <input type="checkbox"/> 15"   | <input type="checkbox"/> 15.5" | <input type="checkbox"/> 16"   | <input type="checkbox"/> 16.5" | <input type="checkbox"/> 17"   | <input type="checkbox"/> 17.5" | <input type="checkbox"/> 18"   | <input type="checkbox"/> 18.5" | <input type="checkbox"/> 19"   | <input type="checkbox"/> 19.5" | <input type="checkbox"/> 20"   | <input type="checkbox"/> 20.5" |                                |                                |  |
|                                       |       | 6"          | <input type="checkbox"/> 15"   | <input type="checkbox"/> 15.5" | <input type="checkbox"/> 16"   | <input type="checkbox"/> 16.5" | <input type="checkbox"/> 17"   | <input type="checkbox"/> 17.5" | <input type="checkbox"/> 18"   | <input type="checkbox"/> 18.5" | <input type="checkbox"/> 19"   | <input type="checkbox"/> 19.5" | <input type="checkbox"/> 20"   | <input type="checkbox"/> 20.5" |                                |                                |                                |  |

### Casters

|             |                          | Caster Size |                          |                          |                          |                          |
|-------------|--------------------------|-------------|--------------------------|--------------------------|--------------------------|--------------------------|
|             |                          | Price       | 3"                       | 4"                       | 5"                       | 6"                       |
| Caster Type | .75" Rollerblade         | N/C         | <input type="checkbox"/> | <input type="checkbox"/> |                          |                          |
|             | .75" Lighted Rollerblade | \$86        | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |                          |
|             | 1.25" Soft Roll          | N/C         |                          | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
|             | 1" Poly Aluminum         | \$131       |                          | <input type="checkbox"/> | <input type="checkbox"/> |                          |
|             | 1" Pneumatic             | \$89        |                          |                          |                          | <input type="checkbox"/> |
|             | 1.5" Soft Roll           | \$68        |                          | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
|             | 1.5" Soft Roll Aluminum  | \$214       | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |

### Fork Options

- ☐ Standard Fork N/C  
☐ Single Sided Fork \$194  
☐ Frog Legs II Suspension System \$508

## Footrest

### Footplate

See reference page (L). Tubular and Flip Under Footrests not available with V-Taper.

- ☐ Tubular "Open" Footplate N/C  
☐ Tubular Footplate with ABS Cover \$32  
☐ Tubular Footplate with Carbon Fiber Cover \$204  
☐ Aluminum Angle Adjustable (K0040) \$102  
☐ Aluminum Angle Adjustable with Riser (K0040) \$172  
☐ Angle Adjustable Flip Under (K0040, K0037) \$236  
☐ Hybrid Angle Adjustable (K0040) \$94  
☐ Hybrid Angle Adjustable with Wings (K0040) \$136

### Footrest Protector

Only available with Tubular Footplates

- ☐ Footrest Protector N/C

### Seat to Footrest Length

See reference page (H). 8" - 9.5" not available with Klaxon Twist.

|                              |                                |                              |                                |                              |                                |
|------------------------------|--------------------------------|------------------------------|--------------------------------|------------------------------|--------------------------------|
| <input type="checkbox"/> 8"  | <input type="checkbox"/> 8.5"  | <input type="checkbox"/> 9"  | <input type="checkbox"/> 9.5"  | <input type="checkbox"/> 10" | <input type="checkbox"/> 10.5" |
| <input type="checkbox"/> 11" | <input type="checkbox"/> 11.5" | <input type="checkbox"/> 12" | <input type="checkbox"/> 12.5" | <input type="checkbox"/> 13" | <input type="checkbox"/> 13.5" |
| <input type="checkbox"/> 14" | <input type="checkbox"/> 14.5" | <input type="checkbox"/> 15" | <input type="checkbox"/> 15.5" | <input type="checkbox"/> 16" | <input type="checkbox"/> 16.5" |
| <input type="checkbox"/> 17" | <input type="checkbox"/> 17.5" | <input type="checkbox"/> 18" | <input type="checkbox"/> 18.5" | <input type="checkbox"/> 19" | <input type="checkbox"/> 19.5" |

## Rear Seat Height and Wheels

### Wheel Size / Rear Seat Height

See reference page (J). Select Rear Seat Height in row corresponding to desired Wheel Size. Selection is required.

| Wheel            | Rear Seat Height               |                                |                                |                                |                                |                                |                                |                                |                                |                                |                              |                                |                              |                                |                              |
|------------------|--------------------------------|--------------------------------|--------------------------------|--------------------------------|--------------------------------|--------------------------------|--------------------------------|--------------------------------|--------------------------------|--------------------------------|------------------------------|--------------------------------|------------------------------|--------------------------------|------------------------------|
| 20" <sup>1</sup> | <input type="checkbox"/> 13"   | <input type="checkbox"/> 13.5" | <input type="checkbox"/> 14"   | <input type="checkbox"/> 14.5" | <input type="checkbox"/> 15"   | <input type="checkbox"/> 15.5" | <input type="checkbox"/> 16"   | <input type="checkbox"/> 16.5" | <input type="checkbox"/> 17"   | <input type="checkbox"/> 17.5" | <input type="checkbox"/> 18" | <input type="checkbox"/> 18.5" | <input type="checkbox"/> 19" | <input type="checkbox"/> 19.5" | <input type="checkbox"/> 20" |
| 22" <sup>2</sup> | <input type="checkbox"/> 14"   | <input type="checkbox"/> 14.5" | <input type="checkbox"/> 15"   | <input type="checkbox"/> 15.5" | <input type="checkbox"/> 16"   | <input type="checkbox"/> 16.5" | <input type="checkbox"/> 17"   | <input type="checkbox"/> 17.5" | <input type="checkbox"/> 18"   | <input type="checkbox"/> 18.5" | <input type="checkbox"/> 19" | <input type="checkbox"/> 19.5" | <input type="checkbox"/> 20" |                                |                              |
| 24" <sup>3</sup> | <input type="checkbox"/> 15"   | <input type="checkbox"/> 15.5" | <input type="checkbox"/> 16"   | <input type="checkbox"/> 16.5" | <input type="checkbox"/> 17"   | <input type="checkbox"/> 17.5" | <input type="checkbox"/> 18"   | <input type="checkbox"/> 18.5" | <input type="checkbox"/> 19"   | <input type="checkbox"/> 19.5" | <input type="checkbox"/> 20" |                                |                              |                                |                              |
| 25" <sup>4</sup> | <input type="checkbox"/> 15.5" | <input type="checkbox"/> 16"   | <input type="checkbox"/> 16.5" | <input type="checkbox"/> 17"   | <input type="checkbox"/> 17.5" | <input type="checkbox"/> 18"   | <input type="checkbox"/> 18.5" | <input type="checkbox"/> 19"   | <input type="checkbox"/> 19.5" | <input type="checkbox"/> 20"   |                              |                                |                              |                                |                              |
| 26" <sup>5</sup> | <input type="checkbox"/> 16"   | <input type="checkbox"/> 16.5" | <input type="checkbox"/> 17"   | <input type="checkbox"/> 17.5" | <input type="checkbox"/> 18"   | <input type="checkbox"/> 18.5" | <input type="checkbox"/> 19"   | <input type="checkbox"/> 19.5" | <input type="checkbox"/> 20"   |                                |                              |                                |                              |                                |                              |

1 Not available with Klaxon Twist

2 Lowest Rear Seat Height with Klaxon Twist and 22" wheels is 15" with Low Profile tires and 15.5" with Full Profile tires

3 Lowest Rear Seat Height with Klaxon Twist and 24" wheels is 16" and 16.5" with Mountain tire

4 Lowest Rear Seat Height with Klaxon Twist and 25" wheels is 16" and 17" with Mountain tire

5 Lowest Rear Seat Height with Klaxon Twist and 26" wheels is 17"

### Optimized Tower

- ☐ Standard (2" Adjustment Range) N/C  
☐ Optimized (1" Adjustment Range) N/C

### Spineray Options

Select spoke color for Spineray Wheels

- ☐ Black Spokes ☐ Red Spokes ☐ Blue Spokes ☐ White Spokes N/C

### Wheels

Select color for Mag (Color) on Color Chart. Superlight wheels available only with Superlight and Plastic Coated Handrims.

|                 |                  | Rear Wheel Size |                          |                          |                          |                          |                          |
|-----------------|------------------|-----------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|
|                 |                  | Price           | 20"                      | 22"                      | 24"                      | 25"                      | 26"                      |
| Rear Wheel Type | MAXX Spoke       | N/C             | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
|                 | Mag              | N/C             | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |                          |                          |
|                 | Mag (Color)      | \$122           | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |                          |                          |
|                 | Superlight Spoke | \$340           |                          |                          | <input type="checkbox"/> | <input type="checkbox"/> |                          |
|                 | Spineray Spox    | \$1014          |                          | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
|                 | Spineray LX      | \$1273          |                          | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |

### Rear Tire Type

Poly, Airless, or Solid tires not available with Superlight wheel. Poly tires not available with 25" or 26" wheels. Schwalbe One and Mountain tires available only with 24" and 25" wheels. Low Profile Schwalbe tires and High Pressure tire not available with Mag wheels. Marathon tire not available with 20" or 22" wheels. Pneumatic tires not available with 25" wheel. Schwalbe HS-127 tire available only with 24" wheels. Low Poly and Sentinel tires not available with Klaxon Twist.

#### Full Profile

- ☐ Street Pneumatic 1 3/8" (E2211, E2212) N/C  
☐ Full Poly 1 3/8" (E2218) N/C  
☐ Pneumatic with Airless Insert 1 3/8" (E2211, E2213) \$101  
☐ Schwalbe HS-127 1 3/8" (E2211, E2212) \$135  
☐ Mountain 2" Black (E2211, E2212) \$153

#### Low Profile

- ☐ Low Poly 1" (E2218) \$81  
☐ Sentinel Solid Tire 1" \$112  
☐ High Pressure 1" (E2211, E2212) \$102  
☐ Schwalbe One 1" (E2211, E2212) \$135  
☐ Schwalbe Marathon Plus Evolution 1" (E2211, E2212) \$167

### Handrims

Superlight handrim not available with 20", 22", or 26" wheels or with Mag wheels. Ergonomic handrims not available with 20" or 26" wheels. Projection handrim not available with Low Profile tires, 20" or 25" Maxx Spoke wheels, or 20" Mag wheels.

#### Standard Handrim

- ☐ Aluminum Anodized N/C  
☐ Black Anodized Superlight \$54  
☐ Plastic Coated \$130  
☐ Projection (8 Vertical Projections) \$374  
☐ Omit Handrim N/C

#### Ergonomic Handrim

- ☐ Ergo Lite Regular Grip \$487  
☐ Ergo Standard Regular Grip \$465  
☐ Ergo Lite Extra Grip \$539  
☐ Ergo Standard Extra Grip \$517  
☐ Ergo Lite Omit Grip \$413  
☐ Ergo Standard Omit Grip \$392

### Handrim Options

#### Handrim Mount

If no selection is made, Standard Mount will be used. Reduces overall width by .75".

- ☐ Close Mount \$28

#### Micro-Grip

Not available on Plastic Coated or Projection Handrims

- ☐ 3M Micro Grip \$68

### Axles

Quad Quick Release not available with Mag or Spineray wheels

- ☐ Quick Release N/C  
☐ Quad Quick Release \$198



Phone (800) 981-1540 | Fax (715) 254-0996  
sales@kimobility.com | www.kimobility.com

# ROGUE DT

## Rear Wheel Spacing

See reference page (K). 1.5" minimum with T-Arms. 1.25" minimum with Spinergy Wheels and Direct/Clamp-On Side Guard Mount.

☐ 0.75" ☐ 1" ☐ 1.25" ☐ 1.5" ☐ 1.75" \_\_\_\_\_ N/C

## Rear Wheel Camber

☐ 0° ☐ 2° ☐ 4° ☐ 6° ☐ 8° \_\_\_\_\_ N/C

## Center of Gravity Preset

See reference page (F)

|                                |                                |                                |                                |                                |
|--------------------------------|--------------------------------|--------------------------------|--------------------------------|--------------------------------|
| <input type="checkbox"/> 1"    | <input type="checkbox"/> 1.25" | <input type="checkbox"/> 1.5"  | <input type="checkbox"/> 1.75" | <input type="checkbox"/> 2"    |
| <input type="checkbox"/> 2.25" | <input type="checkbox"/> 2.5"  | <input type="checkbox"/> 2.75" | <input type="checkbox"/> 3"    | <input type="checkbox"/> 3.25" |
| <input type="checkbox"/> 3.5"  | <input type="checkbox"/> 3.75" | <input type="checkbox"/> 4"    | <input type="checkbox"/> 4.25" | <input type="checkbox"/> 4.5"  |
| <input type="checkbox"/> 4.75" | <input type="checkbox"/> 5"    |                                |                                |                                |

## Wheel Locks

### Wheel Lock

Low Profile not available with Full Poly tires

|                                                   |      |
|---------------------------------------------------|------|
| <input type="checkbox"/> Push to Lock             | N/C  |
| <input type="checkbox"/> Push to Lock Flush Mount | N/C  |
| <input type="checkbox"/> Pull to Lock             | N/C  |
| <input type="checkbox"/> Low Profile Scissor Lock | \$68 |
| <input type="checkbox"/> Omit                     | N/C  |

### Wheel Lock Options

Extension Handles available only with Push to Lock or Pull to Lock Wheel Locks. Grade Aids available only with Push to Lock Wheel Locks and Pneumatic tires.

|                                                    |       |
|----------------------------------------------------|-------|
| <input type="checkbox"/> Extension Handles (E0961) | \$88  |
| <input type="checkbox"/> Grade Aids (E0974)        | \$157 |

## Armrests and Side Guards

### Armrests

Swing Away with Side Guard includes Direct Mount receiver. If T-Arm Low configuration is unavailable, Standard Height T-Arm will be used. Ergonomic Seat Depth not available with T-Arms.

| Armrest Type |                                                | Length                        |                               | Price |
|--------------|------------------------------------------------|-------------------------------|-------------------------------|-------|
|              |                                                | Desk                          | Full                          |       |
|              | Height Adjustable T-Arm Tall (11.5"-15") E0973 | <input type="checkbox"/>      | <input type="checkbox"/>      | \$332 |
|              |                                                | Desk <input type="checkbox"/> | Full <input type="checkbox"/> |       |
|              | Height Adjustable T-Arm (8.5"-13.5") E0973     | <input type="checkbox"/>      | <input type="checkbox"/>      | \$332 |
|              |                                                | Desk <input type="checkbox"/> | Full <input type="checkbox"/> |       |
|              | Height Adjustable T-Arm Low (6.5"-10.5") E0973 | <input type="checkbox"/>      | <input type="checkbox"/>      | \$332 |
|              |                                                | Desk <input type="checkbox"/> | Full <input type="checkbox"/> |       |
|              | Swing Away (8.5"-12.5") E0973                  |                               | <input type="checkbox"/>      | \$249 |
|              | Swing Away with Side Guard (8.5"-12.5") E0973  |                               | <input type="checkbox"/>      | \$298 |
|              | Omit Armrests                                  |                               | <input type="checkbox"/>      | N/C   |

### Side Guard Removal

Selection required with Direct and Clamp-On Side Guard Mount. Not available with Bolt-On and Fixed Side Guard Mount.

|                                        |      |
|----------------------------------------|------|
| <input type="checkbox"/> Removable     | N/C  |
| <input type="checkbox"/> Non-Removable | \$15 |

### Side Guards

Direct Mount with Fender Side Guard not available with Adjustable Depth greater than 1/2". Bolt-On and Fixed Side Guard Mount require Non-Folding Backrest selection.

| Side Guard Type                                                    | Mount Type                                     | Price |
|--------------------------------------------------------------------|------------------------------------------------|-------|
| Composite Side Guards - Adult (K0108)                              | Direct <input type="checkbox"/>                | \$152 |
|                                                                    | Clamp-On <input type="checkbox"/>              | \$183 |
| Aluminum Straight - Black Matte (K0108)                            | Direct <input type="checkbox"/>                | \$204 |
|                                                                    | Clamp-On <input type="checkbox"/>              | \$236 |
| Aluminum Straight - Match Frame Color (K0108)                      | Direct <input type="checkbox"/>                | \$256 |
|                                                                    | Clamp-On <input type="checkbox"/>              | \$288 |
| Aluminum Adjustable Fender w/ Foam Pad - Black Matte (K0108)       | Direct <input type="checkbox"/>                | \$288 |
|                                                                    | Clamp-On <input type="checkbox"/>              | \$319 |
| Aluminum Adjustable Fender w/ Foam Pad - Match Frame Color (K0108) | Direct <input type="checkbox"/>                | \$340 |
|                                                                    | Clamp-On <input type="checkbox"/>              | \$371 |
| Carbon Fiber Fender (K0108)                                        | Direct <input type="checkbox"/>                | \$507 |
|                                                                    | Clamp-On <input type="checkbox"/>              | \$539 |
|                                                                    | Fixed <input type="checkbox"/>                 | \$539 |
| Custom Fit Aluminum Straight - Black Matte (K0108)                 | Bolt-On No Adjustment <input type="checkbox"/> | \$390 |
|                                                                    | Bolt-On 1" Adjustment <input type="checkbox"/> | \$280 |
| Custom Fit Aluminum Straight - Match Frame Color (K0108)           | Bolt-On No Adjustment <input type="checkbox"/> | \$410 |
|                                                                    | Bolt-On 1" Adjustment <input type="checkbox"/> | \$300 |

## Accessories

### Anti-Tips

|                                                 |       |
|-------------------------------------------------|-------|
| <input type="checkbox"/> Rear Anti-Tips (E0971) | \$175 |
| <input type="checkbox"/> User Activated - Pair  | \$356 |
| <input type="checkbox"/> User Activated - Left  | \$177 |
| <input type="checkbox"/> User Activated - Right | \$177 |
| <input type="checkbox"/> Omit                   | N/C   |

### Calf Straps

|                                                            |       |
|------------------------------------------------------------|-------|
| <input type="checkbox"/> Velcro Adjustable (K0038)         | \$44  |
| <input type="checkbox"/> Padded Velcro Adjustable (K0038)  | \$71  |
| <input type="checkbox"/> Padded Airmesh Calf Strap (K0038) | \$123 |

### Pelvic Positioning Belt

|                                                                                    |       |
|------------------------------------------------------------------------------------|-------|
| <input type="checkbox"/> 1 1/2" Auto Buckle Pelvic Positioning Belt (E0978)        | \$91  |
| <input type="checkbox"/> 1 1/2" Padded Auto Buckle Pelvic Positioning Belt (E0978) | \$183 |
| <input type="checkbox"/> Small Flex-Mount Positioning Belt (16-24") (E0978)        | \$208 |
| <input type="checkbox"/> Medium Flex-Mount Positioning Belt (19-29") (E0978)       | \$208 |
| <input type="checkbox"/> Large Flex-Mount Positioning Belt (25-37") (E0978)        | \$208 |

### Other Accessories

|                                                    |       |
|----------------------------------------------------|-------|
| <input type="checkbox"/> Tool Kit                  | \$63  |
| <input type="checkbox"/> Removable Underseat Pouch | \$64  |
| <input type="checkbox"/> Ki Mobility Backpack      | \$83  |
| <input type="checkbox"/> Luggage Carrier           | \$90  |
| <input type="checkbox"/> Cane and Crutch Holder    | \$214 |

### Belt Mounting Kit

*Clamp with bracket defaults when Pelvic Positioning Belt is selected*

|                                                          |       |
|----------------------------------------------------------|-------|
| <input type="checkbox"/> Ki Mobility Belt Mounting Clamp | \$38  |
| <input type="checkbox"/> Belt Mounting Bracket           | \$117 |

## Color

*Color Accent Package includes Backplates (if present), Towers, Tower Clamps, Camber Tube, and Caster Wing Caps. Default accent is black. No selection needed for Backrest Release Bar when Non-Folding Backrest is selected. Touch-up paint available for purchase in the parts manual.*

| Options              | Price | Colors                                  |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |
|----------------------|-------|-----------------------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|
|                      |       | Gloss                                   |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          | Matte                    |                          |                          |                          |                          |
|                      |       | Black                                   | Ruby Red                 | Radiant Orange           | Yellow                   | Electric Green           | Rich Teal                | Glacier Blue             | Blue Steel               | Deep Sapphire Blue       | Purple Plum              | Opalescent Violet        | Dark Pink                | Light Pink               | Champagne                | Shadow Black             | Iridescent Navy          | Bright Red               | Tactical Green           | Dry Earth Tan            |
| Frame                | N/C   | <input type="checkbox"/>                | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| Color Accent Package | N/C   | <input type="checkbox"/>                | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| Backrest Release Bar | N/C   | <input type="checkbox"/> Black Anodized |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |
|                      | \$52  | <input type="checkbox"/>                | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| Mag Color            |       | <input type="checkbox"/>                | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |

### CAD Model

*Visual representation of chair dimensions. If email is not provided, it gets sent to the same email as the quote.*

|                                        |  |
|----------------------------------------|--|
| <input type="checkbox"/> CAD Model     |  |
| <input type="checkbox"/> Email Address |  |

## Power Assist

*The following options are specific to the Twist Power Assist Add-On option*

### Controller Mounting Position

*If no selection is made Right Side Position will be used*

|                                     |     |
|-------------------------------------|-----|
| <input type="checkbox"/> Left Side  | N/C |
| <input type="checkbox"/> Right Side | N/C |

### Installation Options

|                                                               |     |
|---------------------------------------------------------------|-----|
| <input type="checkbox"/> Install Both Front and Rear Hardware | N/C |
| <input type="checkbox"/> Install Front Hardware               | N/C |
| <input type="checkbox"/> Install Rear Hardware                | N/C |
| <input type="checkbox"/> Not Installed                        | N/C |

### Additional Accessories

*Front Hardware required when ordering the Handlebar Controller*

|                                                                        |       |
|------------------------------------------------------------------------|-------|
| <input type="checkbox"/> Attendant Controller                          | \$315 |
| <input type="checkbox"/> Handlebar Controller                          | \$495 |
| <input type="checkbox"/> Twist Carry Bag                               | \$195 |
| <input type="checkbox"/> Handlebar and Front Mount Connector Carry Bag | \$195 |
| <input type="checkbox"/> Removable Li-ion Battery - 36V, 4Ah           | \$595 |

### Hardware Options

*If no selection is made, Rear will be provided. Select only the hardware you need.*

|                                                  |        |
|--------------------------------------------------|--------|
| <input type="checkbox"/> Front and Rear Hardware | N/C    |
| <input type="checkbox"/> Front Hardware          | (\$75) |
| <input type="checkbox"/> Rear Hardware           | (\$75) |

## REFERENCE PAGE

### A. Seat Width (Fig. 1)

Measured from the outside of the horizontal seat tube to the outside of the horizontal seat tube on the opposite side.

### B. Seat Depth (Fig. 2)

Measured from the front of the back post to the front edge of the seat sling.

### C. Sling Position (Fig. 2, Fig. 3)

Measured from the front edge of the seat sling to just before the bend in the seat tube. In the 0", the sling is at the front of the seat tube before the bend. This keeps the chair shorter and more maneuverable. In the 1" or 2" seat sling position the seat tube extends in front of the sling by that indicated amount. Adding more frame can improve stability and provide support to aid in transfers.

### D. Front Frame Bend (Fig. 2)

The angle between the vertical seat tube and the ground.

### E. Ergonomic Seat Depth (Fig. 4)

Measured from the front of the back post forward to where the seat tube begins bending up.

### F. Center of Gravity (Fig. 2)

Measured from the front of the back post to the center of the rear axle.

### G. Footrest Taper (Fig. 1)

Measured from the outside of the horizontal seat tube at the highest point to the outside of the vertical seat tube at the lowest point.

### H. Seat to Footrest Length (Fig. 2)

Measured from the front edge of the seat sling to the footrest. A footrest length of at least 2.5" shorter than the front seat height is recommended.

### I. Front Seat Height (Fig. 2)

Measured from the ground to the top of the seat tube at the front of the seat sling.

### J. Rear Seat Height (Fig. 2)

Measured from the ground to the top of the seat tube at the front of back post.

### K. Rear Wheel Spacing (Fig. 1)

Measured from the outside of the back post to the inside of the rear tire.

### L. Footrest Width (Fig. 1)

Measured from the inside of the vertical seat tube at the lowest point to the inside of the vertical seat tube on the opposite side.

Table 1: Footrest widths are listed below for each seat width and taper.

| Seat Width | Taper    |      |      |       |      |
|------------|----------|------|------|-------|------|
|            | Straight | 1"Y  | 2"Y  | 3.5°V | 7°V  |
| Outside    | 2.25" <  | 4" < | 6" < | 4" <  | 6" < |
| 12         | 9.75     | 8    | N/A  | 8     | N/A  |
| 13         | 10.75    | 9    | N/A  | 9     | N/A  |
| 14         | 11.75    | 10   | 8    | 10    | 8    |
| 15         | 12.75    | 11   | 9    | 11    | 9    |
| 16         | 13.75    | 12   | 10   | 12    | 10   |
| 17         | 14.75    | 13   | 11   | 13    | 11   |
| 18         | 15.75    | 14   | 12   | 14    | 12   |
| 19         | 16.75    | 15   | 13   | 15    | 13   |
| 20         | 17.75    | 16   | 14   | 16    | 14   |




Fig. 1




Fig. 2




Fig. 3




Fig. 4