



Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com

Company: _____

Ship To: _____

Account Number: _____ Date: _____

Address: _____

Quote PO Number: _____

City: _____ State: _____

RTS: _____

Zip: _____ Phone: _____

Marked For: _____

Email/Fax: _____

Select if chair quote has: Additional Ki Products with Different Quote/Order Number **Quote/Order Number:** _____

- Spark - Transit Included** E1238 **\$2395**
 165 lb weight capacity. Not available with upholstery.
- Spark - Transit Not Included** E1238 **\$2075**
 165 lb weight capacity. Not available with upholstery.

Effective 07/18/2022

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions. Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Frame

Frame Width

- 10" N/C
- 11" N/C
- 12" N/C
- 13" N/C
- 14" N/C
- 15" N/C
- 16" N/C

Frame Depth

- 12" N/C
- 13" N/C
- 14" N/C
- 15" N/C
- 16" N/C
- 17" N/C
- 18" N/C

Transit

- Transit System **\$320**
 No charge when selected with Spark - Transit Included. Not available with Removable Seat Pan.
- Omit Transit N/C

Seat Pan

- Removable Aluminum Seat Pan E2231 **\$275**
 Not available with Transit System. Pan depth 1/2" shorter than Frame Depth.
- Omit Seat Pan N/C

Backrest

Back Post / Height Range / Back Height

Back Posts are folding and angle adjustable in 4° increments

Back Post	Price	Height Range	Back Height						
<input type="checkbox"/> Stroller Handle Back Post	N/C	<input type="checkbox"/> Short	<input type="checkbox"/> 20						
		<input type="checkbox"/> Tall	<input type="checkbox"/> 24						
<input type="checkbox"/> Height Adjustable Straight with Push Handle	\$75	<input type="checkbox"/> Short ¹	<input type="checkbox"/> 13	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17		
		<input type="checkbox"/> Medium	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20		
		<input type="checkbox"/> Tall	<input type="checkbox"/> 20 ²	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24		

¹ Not available with Seating Dynamics Rocker Back

² Not available with Angle Adjustable Locking Extendable Flip Up Armrest



Frame Width

Spark is set up to allow 2" increase in frame width from initial frame width in 1" increments. Growth is achieved by expanding cross braces. 16" maximum frame width. 10" Frame Width expands only 1".

Frame Depth

Spark is set up to allow 4" increase in frame depth from initial frame depth in 1" increments. 19" maximum frame depth.

Wheelchair Drive Ready

Chair will ship from Ki Mobility with adaption kit installed for selected Alber system. Order system separately from Alber USA. Select Rear Frame Height (full profile selections only) in drive wheel size on Rear Frame Height chart (page 4). Requires Rear Anti-Tips. Not available with Spinergy wheel.

- e-motion Ready N/C
- twion Ready N/C

Back Post Options

- Removable Stroller Handle Extension **\$75**
 Not available with Height Adjustable Straight with Push Handle Back Posts
- Seating Dynamics Rocker Back E2398 **\$1995**
 Available only with T-Arms. Not available with 12" Frame Depths with Rear Frame Heights lower than 14.5".

Frame Type / Footrest

Frame Type

- Swing Away _____ N/C
- Super Low - Swing Away _____ \$50

Minimum Front Seat Height with Front Mount Hangers						
Hanger Length	Short		Medium		Long	
Hanger Bend	70°	80/90°	70°	80/90°	70°	80°
Minimum Front Seat Height with Standard Selection	13"	13.5"	14"	15"	16"	17"
Minimum Front Seat Height with 3" Extension	14.5"	16"	16.5"	18"	18.5"	20"

Hanger / Front Frame Bend

Hanger	60°	70°	80°	90°		Retail
Extension Tube Mount	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> ¹			Std
Front Mount		<input type="checkbox"/>	<input type="checkbox"/> ¹	<input type="checkbox"/> ¹		N/C
Pro Elevating Leg Rests ² (PRO ELR) E0990					<input type="checkbox"/>	\$295
Pediatric Pro Elevating Leg Rests (Peds PRO ELR) E0990					<input type="checkbox"/>	\$295
Omit Hangers					<input type="checkbox"/>	-\$20

¹ Not available with 7" or 8" Casters

² Not available with Seat Widths less than 14"

PRO ELR Covers

Available only on PRO ELR

- Gel Pad Covers (Pair) K0108 _____ \$390

Footrest / Extension / Seat to Footrest Length

Select Extension and Seat to Footrest Length in row corresponding to desired Footrest. Footrest Length of at least 2.5" shorter than Front Seat Height recommended. Extension Tube Hanger measurements are based on 70° Hanger.

Footrest	Extension	Seat to Footrest Length											
Extension Tube Hangers (60°, 70°, 80°)	<input type="checkbox"/> Short	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15			
	<input type="checkbox"/> Medium	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17			
	<input type="checkbox"/> Long	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5			
Front Mount Hangers ¹	<input type="checkbox"/> Short (70°, 80°, 90°)	Standard Selections									Requires 3" Extension for \$55		
		<input type="checkbox"/> 4 ⁵	<input type="checkbox"/> 4.5 ⁵	<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5		<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11		
	<input type="checkbox"/> Medium (70°, 80°, 90°)	<input type="checkbox"/> 7	<input type="checkbox"/> 7.5	<input type="checkbox"/> 8	<input type="checkbox"/> 8.5	<input type="checkbox"/> 9	<input type="checkbox"/> 9.5		<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5		
		<input type="checkbox"/> 4 ⁵	<input type="checkbox"/> 4.5 ⁵	<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5	<input type="checkbox"/> 7	<input type="checkbox"/> 7.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	
	<input type="checkbox"/> Long (70°, 80°)	<input type="checkbox"/> 8	<input type="checkbox"/> 8.5	<input type="checkbox"/> 9	<input type="checkbox"/> 9.5	<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	
		<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5	<input type="checkbox"/> 7	<input type="checkbox"/> 7.5	<input type="checkbox"/> 8	<input type="checkbox"/> 8.5	<input type="checkbox"/> 9	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15
Peds ² PRO ELR	<input type="checkbox"/> Low Pivot	<input type="checkbox"/> 7.5 ³	<input type="checkbox"/> 8 ³	<input type="checkbox"/> 8.5 ³	<input type="checkbox"/> 9 ⁴	<input type="checkbox"/> 9.5 ⁴	<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13
	<input type="checkbox"/> High Pivot	<input type="checkbox"/> 7.5 ³	<input type="checkbox"/> 8 ³	<input type="checkbox"/> 8.5 ³	<input type="checkbox"/> 9 ⁴	<input type="checkbox"/> 9.5 ⁴	<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5			
PRO ELR Low Pivot Setting	<input type="checkbox"/> Short	12 - 17.5											
	<input type="checkbox"/> Long	15 - 21											
PRO ELR High Pivot Setting	<input type="checkbox"/> Short	10.5 - 15.5											
	<input type="checkbox"/> Long	13 - 19											

¹ See Minimum Front Seat Height chart above.

² Comes standard with 3" Extension

³ Not available on chairs 13" wide or more

⁴ Not available on chairs 15" wide or more

⁵ Available with 90° only

3" Extension

Available only with Front Mount Hangers

- 3" Extension _____ \$55
- 3" Extension (Provided as Parts) _____ \$55


Footplate

- Composite _____ N/C
Not available with Front Mount Hanger or Frame Widths less than 14"
- Composite Angle Adjustable K0040 _____ \$190
Not available with Front Mount Hanger or Frame Widths less than 15"
- Aluminum Angle Adjustable K0040 _____ \$190
- Aluminum Locking Multi-Angle Adjustable K0108 _____ \$300
Available only with Front Mount Hanger
- One Piece Flip Up Angle Adjustable K0040 _____ \$190
Not available with Front Mount Hanger, PRO ELR or Frame Widths less than 14"

Hanger Release

If no selection is made Classic Hanger Release will be used

- Classic _____ Std
- 4-Way - In line _____ N/C
See reference A below. Not available with PRO ELR or Peds PRO ELR.
- 4-Way - Out _____ N/C
See reference B below. Not available with PRO ELR or Peds PRO ELR.



Frame Type / Footrest continued

Heel Loop

Not available with One Piece Flip Up Angle Adjustable Footplate

- Heel Loops E0951 _____ \$35
Not available with Frame Widths less than 14"
- Adjustable Heel Loops E0951 _____ \$50
- Omit Heel Loops _____ N/C

Calf Strap

- Velcro Adjustable K0038 _____ \$40
- Padded Velcro Adjustable K0038 _____ \$68
- Bodypoint Padded AeroMesh K0038 _____ \$118

Residual Limb Support

Residual Limb Support E1020

- Right and Left _____ \$810
- Left _____ \$405 Right _____ \$4055
- Omit Left Hanger _ -\$10 Omit Right Hanger _ -\$10

Cover

- Smooth Air Mesh Smooth Air Mesh

Pad Size

- | | | | |
|---------------------------------|----------------------------------|---------------------------------|----------------------------------|
| <input type="checkbox"/> 5 x 8 | <input type="checkbox"/> 9 x 14 | <input type="checkbox"/> 5 x 8 | <input type="checkbox"/> 9 x 14 |
| <input type="checkbox"/> 7 x 10 | <input type="checkbox"/> 11 x 10 | <input type="checkbox"/> 7 x 10 | <input type="checkbox"/> 11 x 10 |
| <input type="checkbox"/> 9 x 10 | <input type="checkbox"/> 11 x 14 | <input type="checkbox"/> 9 x 10 | <input type="checkbox"/> 11 x 14 |

Front Frame Height and Casters

Caster Size / Front Frame Type / Front Frame Height

Select Front Frame Height in row corresponding to desired caster size. Selections in first three columns require Super Low Front Frame.

Caster Size	Front Frame Height													
	Requires Super Low Front Frame ¹													
3"	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14		<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5		
4"	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20
5"	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	
6"	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5			
7"	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5		<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5					
8"	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5							

¹ Price is \$50

Caster

Caster Type		Caster Size						
		Price	3"	4"	5"	6"	7"	8"
Caster Type	.75" Rollerblade	N/C	<input type="checkbox"/>	<input type="checkbox"/>				
	.75" Lighted Rollerblade	\$70	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
	1" Poly	N/C		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	1" Poly Aluminum	\$110		<input type="checkbox"/>	<input type="checkbox"/>			
	1" Pneumatic	\$90				<input type="checkbox"/>		<input type="checkbox"/>
	1.5" Poly	\$50		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	1.5" Soft Roll Aluminum	\$185	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	2" Poly	\$70				<input type="checkbox"/>		
2" Pneumatic w/ Foam Insert	\$225						<input type="checkbox"/>	

Fork Stem

- Best Fit _____ Std
Ki Mobility will pick Fork and Stem combination to optimize performance based on Caster Size and Front Seat Height
- Standard Fork Stem _____ N/C
- 3/4" Longer Fork Stem _____ N/C
- 1 1/2" Longer Fork Stem _____ N/C

Maximum Seat Dump Reference

(Front Seat Height minus Rear Seat Height)

Seat Depth	12"	13"	14"	15"	16"	17"	18"	19"
Maximum Dump Standard Caster Housing	1.5"	1.5"	1.5"	2"	2"	2"	2.5"	2.5"
Maximum Dump Performance Caster Housing	2.5"	2.5"	3"	3"	3"	3.5"	3.5"	4"

Caster Housing

- Standard _____ N/C
- Performance _____ \$50
Not available with 350 lb Weight Capacity

Fork Options

- Frog Legs Suspension Fork _____ \$425
Available only with 3", 4", 5" or 6" Caster Sizes. Not available with Caster Pin Locks, 2" Poly Caster Type or 2" Pneumatic w/ Foam Insert Caster Type.
 - User Weight Under 80 lbs
 - User Weight 80 - 120 lbs
 - User Weight 120 - 160 lbs
- Caster Pin Locks K0073 _____ \$195
Not available with Performance Caster Housing, 1 1/2" Longer Fork Stem or Frog Legs Fork

Wheel Locks

Wheel Lock

Push to Lock _____ Std

Pull to Lock _____ N/C

Omit Wheel Locks _____ N/C

Wheel Lock Options

Extension Handles E0961 _____ \$75
Available only with Push to Lock or Pull to Lock Wheel Locks

Grade Aids E0974 _____ \$150
Available only with Push to Lock Wheel Locks and Pneumatic and Pneumatic with Airless Insert tires

Rear Frame Height and Wheels

Wheel Size / Tire Profile / Rear Frame Height

Select Rear Frame Height in row corresponding to Wheel Size and Tire Profile. 12"-14" Frame Depths not available with 24" wheel size.

Requires Universal Axle Plate Any axle plate can be used (Universal, Amputee, One Arm Drive)

Wheel Size	Tire Profile	Rear Frame Height										
		<input type="radio"/> 11	<input type="checkbox"/> 11.5 ¹	<input type="radio"/> 12	<input type="checkbox"/> 12.5 ¹	<input type="radio"/> 13	<input type="checkbox"/> 13.5	<input type="radio"/> 14	<input type="checkbox"/> 14.5	<input type="radio"/> 15	<input type="radio"/> 15.5	<input type="radio"/> 16
18"	Low	<input type="radio"/> 11	<input type="checkbox"/> 11.5 ¹	<input type="radio"/> 12	<input type="checkbox"/> 12.5 ¹	<input type="radio"/> 13	<input type="checkbox"/> 13.5	<input type="radio"/> 14	<input type="checkbox"/> 14.5	<input type="radio"/> 15	<input type="radio"/> 15.5	<input type="radio"/> 16
	Full	<input type="radio"/> 11.5	<input type="checkbox"/> 12 ¹	<input type="radio"/> 12.5	<input type="checkbox"/> 13 ¹	<input type="radio"/> 13.5	<input type="checkbox"/> 14	<input type="radio"/> 14.5	<input type="checkbox"/> 15	<input type="radio"/> 15.5	<input type="radio"/> 16	<input type="radio"/> 16.5
20"	Low	<input type="radio"/> 11.5	<input type="checkbox"/> 12 ¹	<input type="radio"/> 12.5	<input type="checkbox"/> 13 ¹	<input type="radio"/> 13.5	<input type="checkbox"/> 14	<input type="radio"/> 14.5	<input type="checkbox"/> 15	<input type="radio"/> 15.5	<input type="radio"/> 16	<input type="radio"/> 16.5
	Full	<input type="radio"/> 12	<input type="checkbox"/> 12.5 ¹	<input type="radio"/> 13	<input type="checkbox"/> 13.5 ¹	<input type="radio"/> 14	<input type="checkbox"/> 14.5	<input type="radio"/> 15	<input type="checkbox"/> 15.5	<input type="radio"/> 16	<input type="radio"/> 16.5	<input type="radio"/> 17
22"	Low	<input type="radio"/> 12.5	<input type="checkbox"/> 13 ¹	<input type="radio"/> 13.5	<input type="checkbox"/> 14 ¹	<input type="radio"/> 14.5	<input type="checkbox"/> 15	<input type="radio"/> 15.5	<input type="checkbox"/> 16	<input type="radio"/> 16.5	<input type="radio"/> 17	<input type="radio"/> 17.5
	Full	<input type="radio"/> 13	<input type="checkbox"/> 13.5 ¹	<input type="radio"/> 14	<input type="checkbox"/> 14.5 ¹	<input type="radio"/> 15	<input type="checkbox"/> 15.5	<input type="radio"/> 16	<input type="checkbox"/> 16.5	<input type="radio"/> 17	<input type="radio"/> 17.5	<input type="radio"/> 18
24"	Low	<input type="radio"/> 13.5	<input type="checkbox"/> 14 ¹	<input type="radio"/> 14.5	<input type="checkbox"/> 15 ¹	<input type="radio"/> 15.5	<input type="checkbox"/> 16	<input type="radio"/> 16.5	<input type="checkbox"/> 17	<input type="radio"/> 17.5	<input type="radio"/> 18	<input type="radio"/> 18.5
	Full	<input type="radio"/> 14	<input type="checkbox"/> 14.5 ¹	<input type="radio"/> 15	<input type="checkbox"/> 15.5 ¹	<input type="radio"/> 16	<input type="checkbox"/> 16.5	<input type="radio"/> 17	<input type="checkbox"/> 17.5	<input type="radio"/> 18	<input type="radio"/> 18.5	<input type="radio"/> 19

¹ Amputee Axle Plate not available with lowest two Rear Frame Height positions

One Arm Drive

One Arm Drive E0958 _____ \$1095
*Available only with 20", 22" and 24" spoke wheels and Quick Release Axles.
 Spoke wheel style is exclusive to OAD. Not available with 4 degrees of camber.
 Standard handrim will be Aluminum Anodized.*

Rear Wheel Size: 20" 22" 24"

Drive Handrim Side: Left Right

Drive Handrim Type: Aluminum _____ N/C
 Poly Coated _____ \$200

Axle Plate

No selection needed with One Arm Drive

Universal _____ Std

Amputee _____ \$175
Not available with Wheelchair Drive Ready Option

Rear Wheel

No selection needed with One Arm Drive

Rear Wheel Type	MAXX Performance Spoke	Rear Wheel Size				
		Price	18"	20"	22"	24"
	<input type="checkbox"/> MAXX Performance Spoke	\$100	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	<input type="checkbox"/> Mag (Black)	N/C		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	<input type="checkbox"/> Colored Mag	\$100		<input type="checkbox"/>	<input type="checkbox"/>	
	<input type="checkbox"/> Spinergy Spox	\$845			<input type="checkbox"/>	<input type="checkbox"/>

Spinergy Options

Select spoke color and handrim tab for Spinergy wheels

Spoke Color: Black Red White Blue Yellow

Colored Mag

Mag Color: Florescent Red White Blue Pink

Handrim

No selection needed with One Arm Drive

Standard Handrims

Aluminum Anodized _____ Std

Plastic Coated _____ \$95

Projection (8 Vertical Projections) _____ \$325
Available only with full profile tires in combination with 22"-24" MAXX Performance or 24" Mag wheel

Omit Handrim _____ N/C

Natural Fit Handrims

Not available on 18" or 20" wheel sizes or 22" Mag wheel. Not recommended with Mag wheel.

Natural Fit LT Regular Grip _____ \$465

Natural Fit Standard Regular Grip _____ \$445

Natural Fit LT Super Grip _____ \$515

Natural Fit Standard Super Grip _____ \$495

Natural Fit LT Omit Grip _____ \$395

Natural Fit Standard Omit Grip _____ \$345

Handrim Options

3M Micro Grip _____ \$65
Not available on Plastic Coated or Projection Handrims

Handrim Tab:

Available only with Spinergy, Mag or Maxx Performance Spoke and Aluminum Anodized or Plastic Coated Handrims

Short Tab _____ N/C

Long Tab _____ N/C

Rear Tire Type

Full Profile

Street Pneumatic 1 3/8" E2211, E2212 _____ Std

Full Poly 1 3/8" E2218 _____ \$55

Pneumatic with Airless Insert 1 3/8" E2211, E2213 _____ \$95

Low Profile

Low Poly 1" E2218 _____ \$55

Sentinel Solid Tire 1" E2218 _____ \$95
Not available on 18" wheel size or 20" and 22" Mag wheels

High Pressure 1" E2211, E2212 _____ \$95
Not available on Mag wheel or 18" wheel size

Axles


Axle	Axle Position
No selection needed with One Arm Drive	
<input type="checkbox"/> Quick Release _____ Std	Rear Wheel Camber: <input type="checkbox"/> 0° <input type="checkbox"/> 2° <input type="checkbox"/> 4°
<input type="checkbox"/> Quad Quick Release _____ \$175 <i>Not available with Mag or Spinerigy wheels</i>	Center of Gravity: <input type="checkbox"/> Rear <input type="checkbox"/> Center <input type="checkbox"/> Forward <i>No selection needed with Amputee Axle Plate</i>

Armrests

Armrest	Arm Pad
Height Adjustable T-Arm Tall (11.5" - 15") E0973 New!	Angle Adjustable Locking Extendable Flip Up E0973
<input type="checkbox"/> Desk _____ \$310	<input type="checkbox"/> Desk _____ \$265
<input type="checkbox"/> Full _____ \$310	<input type="checkbox"/> Full _____ \$265
Height Adjustable T-Arm (8.5" - 13.5") E0973	<input type="checkbox"/> Pediatric Height Adjustable T-Arm Desk (6.5" - 10") E0973 _____ \$195
<i>Not available with Frame Depths less than 14"</i>	<input type="checkbox"/> Tubular Flip Up _____ N/C <i>Not available with Arm Pads</i>
<input type="checkbox"/> Desk _____ \$310	<input type="checkbox"/> Omit Armrests _____ N/C
<input type="checkbox"/> Full _____ \$310	
Height Adjustable T-Arm Low (6.5" - 10.5") E0973	
<i>Not available with Frame Depths less than 14". If requested configuration is unavailable, Standard Height T-Arm will be used.</i>	
<input type="checkbox"/> Desk _____ \$310	
<input type="checkbox"/> Full _____ \$310	

Accessories

Removable Side Guard	Other Accessories
Not available with Height Adjustable T-Arm Armrests	
<input type="checkbox"/> Adult Composite Side Guards (14" - 16") K0108 _____ \$145	<input type="checkbox"/> Tool Kit _____ \$45
<input type="checkbox"/> Pediatric Composite Side Guards (10" - 14") K0108 _____ \$145	<input type="checkbox"/> Removable Underseat Pouch _____ \$50 <i>Not available with 10" - 13" Seat Widths</i>
Anti-Tip	<input type="checkbox"/> Pediatric Ki Mobility Backpack _____ \$65
<input type="checkbox"/> Rear Anti-Tips E0971 _____ \$145	<input type="checkbox"/> Ki Mobility Backpack _____ \$65
<input type="checkbox"/> Omit Rear Anti-Tips _____ N/C	<input type="checkbox"/> Black Neoprene Impact Guards _____ \$90
Pelvic Positioning Belt	<input type="checkbox"/> Cane and Crutch Holder _____ \$195
<input type="checkbox"/> 1" Auto Buckle E0978 _____ \$80	<input type="checkbox"/> Canopy with Hardware K0108 _____ \$280 <i>See image below</i>
<input type="checkbox"/> 1" Padded Auto Buckle E0978 _____ \$170	
<input type="checkbox"/> 1 1/2" Auto Buckle E0978 _____ \$80	
<input type="checkbox"/> 1 1/2" Padded Auto Buckle E0978 _____ \$170	
<input type="checkbox"/> 2" Airline Buckle E0978 _____ \$110	
Bodypoint Evoflex - Pelvic Stabilizer E0978 _____ \$188	
Size: <input type="checkbox"/> Small (16-24") <input type="checkbox"/> Medium (19-29") <input type="checkbox"/> Large (25-37")	
Bodypoint Belt Mounting Kit _____ \$108	
Style: <input type="checkbox"/> Bracket <input type="checkbox"/> Band Clamp 1"	
Spoke Guards	
<input type="checkbox"/> Clear Spoke Guards K0065 _____ \$115 <i>Not available with Mag wheel</i>	
<input type="checkbox"/> Black Spoke Guards K0065 _____ \$115 <i>Not available with 18" or 20" wheel sizes or Mag wheel</i>	
Designer Spoke Guards K0065 _____ \$250 <i>Available only with 20" and 22" wheel sizes. Not available with Mag wheel.</i>	
Design:	
<input type="checkbox"/> Blue Monster <input type="checkbox"/> Fairy Garden <input type="checkbox"/> Havin' a Hoot	
<input type="checkbox"/> Robo Fun <input type="checkbox"/> Rockin' Reptiles <input type="checkbox"/> Super Puppies	
<input type="checkbox"/> Tutu Cute <input type="checkbox"/> Tuxedo Time	



Frame Color

Gloss Finish

- Black _____ N/C
- Blue Skies _____ N/C
- Bubblegum Pink _____ N/C
- Candy Blue _____ N/C
- Candy Purple _____ N/C
- Candy Red _____ N/C
- Copper _____ N/C
- Granny Smith _____ N/C
- Metallic Orange _____ N/C
- Yellow _____ N/C

Matte Finish

- Charcoal Gray _____ N/C
- Emerald Green _____ N/C
- Iridescent Burgundy _____ N/C
- Iridescent Navy _____ N/C
- Military Green _____ N/C
- Retro Red _____ N/C
- Royal Blue _____ N/C
- Shadow Black _____ N/C
- Smooth White _____ N/C

Hydrographics

Hydrographics _____ \$775

- AOR1 (Desert Digital Camo)
- Pink Camo
- Muddy Girl
- NEXT Camo
- Carbon Fiber¹
- Zebra¹

¹Select Hydrographics Base Color. Zebra not available in Black.

Hydrographics Base Color _____ N/C

- Black
- Blue
- Granny Smith
- Pink
- Purple
- Red

Spark Front Frame Height Matrix																																			
3" Caster				4" Caster				5" Caster				6" Caster				7" Caster				8" Caster															
FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork												
14.5	13.5	Std	4"	15	13.5	Std	4"	16	14.5	Std	4"	17	15.5	Std	5"	18	16.5	Std	6"	19	17.5	Std	6"												
15	14			15.5	14			16.5	15.5			3/4	17.5			16.5	3/4			18.5	17			19.5	18.5	3/4									
15.5	14	3/4		15.5 ¹	14 ¹	3/4		17.5	16	1 1/2		18.5	17	1 1/2		18.5	17.5	3/4		20.5	n/a	1 1/2													
16	14.5			16	14.5			Std	17	15.5		1 1/2	17.5	16		Std	19			18	1 1/2	19		17.5	Std										
16	15	1 1/2	16.5 ¹	15 ¹	1 1/2	16.5	15		3/4	17.5	16		5"	17.5	16		Std	19.5	18	Std															
16.5	15.5		17	15.5		3/4	16.5	15.5		5"	18	16.5		6"	18	16.5		Std	20		18.5	Std													
15.5	14.5	Std	15.5	14	Std		17	16	1 1/2		17.5	16	3/4		17.5	16.5	6"		18.5	17	Std														
16	15		3/4	16		14.5	3/4	17.5		16	5"	18		17	Std	18		17	Std	19		17.5	Std												
16.5	15	1 1/2		16	15	3/4		18	16.5	Std		18.5	17.5	6"		18.5	17.5	Std		19.5	18	Std													
17	15.5		1 1/2	16.5	15.5		1 1/2	16.5	15		Std	18.5	17		7"	18.5	17		3/4	20	18.5		Std												
17	16	Std		17	15.5	Std		17	15.5	3/4		19	17.5	6"		19	17.5	Std		19	18	Std													
17.5	16.5		3/4	17.5	16		3/4	17.5	16		1 1/2	19.5	18		7"	19.5	18		Std	19.5	18.5		Std												
16.5	15.5	1 1/2		16.5	15	1 1/2		17.5	16	Std		18	16.5	3/4		18	16.5	Std		20	19	Std													
17	16		Std	17	15.5		Std	18	16.5		3/4	18.5	17		6"	18.5	17		Std	20.5	n/a		Std												
17.5	16	3/4		17	16	3/4		18.5	17	1 1/2		19	17.5	7"		19	17.5	Std		20	18.5	Std													
18	16.5		1 1/2	17.5	16.5		1 1/2	18	16.5		Std	18.5	17		3/4	19.5	18		1 1/2	20.5	19		Std												
18	17	Std		18	16.5	Std		18.5	17	3/4		18.5	17	6"		18.5	17.5	Std		21	n/a	Std													
18.5	17.5		1 1/2	18.5	17		1 1/2	19	17.5		Std	19	18		7"	19	18		3/4	20.5	19		Std												
			18				16.5				Std				18				16.5				Std						19.5	18.5	3/4	19.5	18.5	7"	6" x 2" caster not available with 5" fork
			18.5				17								18.5				17										20	19					
			18.5				17.5				3/4				19				17.5				3/4						19.5	18	1 1/2				
			19				18								18.5				17.5										20	18.5					
			19.5				18.5				1 1/2				19				18				3/4						20.5	19	7"				
20	19	19.5	18.5	21	n/a																														
			19.5	18	1 1/2	20	18.5	7"																											
			20	18.5		20.5	19																												
			20.5	19																															

¹ Not available with 4 x 1.5" Soft Roll Aluminum

Seat heights with "FSH" heading available with Standard Front Frame. Seat heights with "Super Low" heading available with Super Low Front Frame. Actual seat height may vary depending upon configuration.