



Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com

TSUNAMI AL

Company: _____

Ship To: _____

Account Number: _____ Date: _____

Address: _____

Quote PO Number: _____

City: _____ State: _____

RTS: _____

Zip: _____ Phone: _____

Marked For: _____

Email/Fax: _____

Select if chair quote has: Additional Ki Products with Different Quote/Order Number **Quote/Order Number:** _____

Effective 01/17/2022

Tsunami AL K0005 \$2455
 250 lb weight capacity

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions. Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Frame

Seat Width

See reference A on page 6

- 14" _____ N/C
- 15" _____ N/C
- 16" _____ N/C
- 17" _____ N/C
- 18" _____ N/C
- 19" _____ N/C
- 20" E2201 \$500
- 21" E2201 \$500
- 22" E2201 \$500

Seat Depth

See reference B on page 6

- 14" _____ N/C
- 15" _____ N/C
- 16" _____ N/C
- 17" _____ N/C
- 18" _____ N/C
- 19" _____ N/C
- 20" E2203 \$500

Weight Capacity

- 250 lb _____ N/C
- 350 lb (Heavy Duty) _____ \$750
Requires 6" Casters

Transit

- Transit System _____ \$320
125 lb maximum weight capacity
- Omit Transit _____ N/C

Backrest

Back Post / Height Range / Back Height

Back Post	Price	Height Range	Back Height				
Height Adjustable	N/C	<input type="checkbox"/> Short	<input type="checkbox"/> 11 ²	<input type="checkbox"/> 12	<input type="checkbox"/> 13	<input type="checkbox"/> 14	
		<input type="checkbox"/> Medium	<input type="checkbox"/> 13	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17
		<input type="checkbox"/> Tall	<input type="checkbox"/> 17	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20	
Fixed Height ¹	\$50	<input type="checkbox"/> Fixed	<input type="checkbox"/> 9 ²	<input type="checkbox"/> 10 ²	<input type="checkbox"/> 11 ²	<input type="checkbox"/> 12 ³	<input type="checkbox"/> 13
			<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18

¹ Not available with back upholstery

² Not available with Fold Down Push Handle

³ Fold Down Push Handle requires Height Adjustable Rigidizer Bar

Sling Position

In Performance Position (0"), sling is at front of frame. The 1" or 2" sling position extends the frame in front of the sling by that amount. See reference C on page 6.

- 0" (Performance) _____ N/C
- 1" _____ N/C
- 2" (Classic) _____ N/C

Adjustable Depth

Adds the capability to grow seat in depth

- No Rear Adjustment _____ N/C
- 1/2" _____ N/C
- 1" _____ N/C
- 2" _____ \$165
- 3" _____ \$165

Long Sling

Available only with 1", 2" and 3" Adjustable Depth

- Long Sling _____ N/C
Sling will be maximum 2" longer than seat depth

Seat Option

- Standard Seat Sling _____ N/C
- Solid Seat Pan E2231 _____ \$265
Seat Pan is 1" longer when Long Sling option is selected

Front Frame Bend

Angle between front seat tube and ground. See reference D on page 6.

- 90° _____ \$75
- 85° _____ N/C
- 80° _____ N/C
- 75° _____ N/C
- 70° _____ N/C

Push Handle

- Omit Push Handle _____ N/C
- Standard Push Handle _____ \$85
- Bolt-On Push Handle _____ \$125
Handles will be 1" - 2" above back height
- Fold Down Push Handle _____ \$100
Not available with 350 lb Weight Capacity

Backrest continued

Stroller Handle

- Stroller Handles are removable and height adjustable*
- Single Stroller Handle _____ \$165
 - Double Stroller Handle _____ \$330
Not available with Push Handles

Rigidizer Bar

- Adjustable Height Rigidizer Bar _____ \$75
- Non-Adjustable Height Rigidizer Bar _____ N/C

Seat Back Angle

Standard is 90° (perpendicular to floor). Minus (-°) closes back angle. Plus (+°) opens back angle.

- Standard +3° +6° +9°
- 3° -6° -9° -12°

Back Upholstery

Available only with Height Adjustable Back Posts. If selecting both upholstery and Axiom Back, the upholstery will be installed.

- Padded Velcro Adjustable
Not available with 350 lb Weight Capacity
 - Standard Cover _____ Std
 - Air Mesh Breathable Cover _____ \$75
- Tension Adjustable E2611
 - Standard Cover _____ \$300
 - Air Mesh Breathable Cover _____ \$375
- Omit Back Upholstery _____ - \$10

Military Patch

- US Military Service Patch _____ \$90
 - Army Navy Marines Air Force Coast Guard
- US Purple Heart Patch _____ \$90
- Purple Heart with Military Service Patch Combo _____ \$145
 - Army Navy Marines Air Force Coast Guard

Front Seat Height and Casters

Caster Size / Footrest Taper / Front Seat Height

Select Front Seat Height and Footrest Taper in row corresponding to Caster Size. 6" casters not available with 85° and 90° Front Frame Bend. Footrest Taper indicates bend of front frame creating footrest width. Taper is measured from outside of seat frame to outside of front tube. 2" Taper not available with 75° Front Frame Bend. Taper is not available with 70° Front Frame Bend. See reference E on page 6. Front Seat Height is measured from floor to top of seat tube at front of seat sling. Adjustable up or down by 1/2". See reference F on page 6.

Footrest Taper	Caster Size	Front Seat Height														
<input type="checkbox"/> Straight	3"	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5
	4"	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	
	5"	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5		
	6"	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5			
<input type="checkbox"/> 1" Taper	3"	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5						
	4"	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5							
	5"	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5								
	6"	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5									
<input type="checkbox"/> 2" Taper	3"	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5						
	4"	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5							
	5"	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5								
	6"	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5									

Caster

Caster Type		Caster Size				
		Price	3"	4"	5"	6"
Caster Type	.75" Rollerblade	N/C	<input type="checkbox"/>	<input type="checkbox"/>		
	.75" Lighted Rollerblade	\$70	<input type="checkbox"/>	<input type="checkbox"/>		
	1" Poly	N/C		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	1" Poly Aluminum	\$110		<input type="checkbox"/>	<input type="checkbox"/>	
	1" Pneumatic	\$80				<input type="checkbox"/>
	1.5" Poly	\$55		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	1.5" Soft Roll Aluminum	\$185	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Fork Options

- Slotted Fork _____ Std
- Single Sided Fork _____ \$170
Not available with 6" Casters
- Frog Legs Suspension Fork _____ \$450
Not available with 350lb Weight Capacity
 - User Weight 80 - 120 lbs
 - User Weight 120 - 160 lbs
 - User Weight 160 - 220 lbs
 - User Weight 220 - 250 lbs
- Frog Legs II Suspension System _____ \$450
Lightweight shock absorbing fork. Not available with 350 lb Weight Capacity.
- Caster Pin Locks K0073 _____ \$185
Available only with Slotted Fork

Footrest

- Footplate**
- Tubular "Open" Footplate _____ Std
Not available w/footrest widths greater than 17.5". See Reference G on page 6.
 - Tubular Footplate with ABS Cover _____ \$30
Not available w/footrest widths greater than 15.5". See Reference G on page 6.
 - Aluminum Angle Adjustable K0040 _____ \$95
 - Aluminum Angle Adjustable with Riser K0040, K0037 _____ \$145
Not available with 350 lb Weight Capacity
 - High Mount Aluminum Angle Adjustable K0040, K0037 _____ \$95
Not available with 350 lb Weight Capacity
 - Angle Adjustable Flip Under K0040, K0037 _____ \$185
Not available with footrest widths less than 9.5". Not available with 350 lb Weight Capacity.
 - High Mount Angle Adjustable Flip Under K0040, K0037 _____ \$270
Not available with footrest widths greater than 15.5" or 350 lb Weight Capacity. See Reference G on page 6.

- Footrest Protector**
- Footrest Protector _____ \$25
Not available with Angle Adjustable Flip Under, High Mount Angle Adjustable Flip Under or High Mount Aluminum Angle Adjustable Footplates

Seat to Footrest Length

Measured from front edge of seat sling to top, rear of footrest. Footrest length at least 2.5" shorter than Front Seat Height recommended. See reference H on page 6.

<input type="checkbox"/> 8"	<input type="checkbox"/> 8.5"	<input type="checkbox"/> 9"	<input type="checkbox"/> 9.5"	<input type="checkbox"/> 10"	<input type="checkbox"/> 10.5"
<input type="checkbox"/> 11"	<input type="checkbox"/> 11.5"	<input type="checkbox"/> 12"	<input type="checkbox"/> 12.5"	<input type="checkbox"/> 13"	<input type="checkbox"/> 13.5"
<input type="checkbox"/> 14"	<input type="checkbox"/> 14.5"	<input type="checkbox"/> 15"	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"
<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"

Rear Seat Height and Wheels

Wheel Size / Rear Seat Height

Select Rear Seat Height in row corresponding to desired Wheel Size. For proper setup, selection is required when omitting rear wheels. Measured from floor to top of seat tube at front of back post. Rear Seat Height is custom to needs. See reference I on page 6.

Wheel Size	Rear Seat Height															
20"	<input type="checkbox"/> 12.5 ¹	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20
22"	<input type="checkbox"/> 13.5 ¹	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20		
24"	<input type="checkbox"/> 14.5 ¹	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20				
25"	<input type="checkbox"/> 15 ¹	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20					
26"	<input type="checkbox"/> 15.5 ¹	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20						

¹Available only with Low Profile tires

- Camber Bar Tower**
- Standard Tower _____ N/C
 - Mx Tower _____ \$100
- Standard Tower** allows rear seat height adjustment up or down in 1/4" increments. 1" adjustment available; 1/2" in lowest end of range. **Mx Tower** required with select low seat heights or when large rear seat height adjustment range is required. Provides 3" adjustment range in 1/4" increments.

Rear Wheel

Rear Wheel	Price	Rear Wheel Size				
		20"	22"	24"	25"	26"
MAXX Performance Spoke	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Mag	\$50	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
Superlight Spoke ¹	\$295			<input type="checkbox"/>	<input type="checkbox"/> ²	
Spinergy Spox	\$845		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Spinergy LX	\$1050		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

¹ Not available with 350 lb Weight Capacity
² Available only with Superlight and Plastic Coated Handrims

- Spinergy Options**
- Select spoke color and handrim tab for Spinergy wheels*
- Spoke Color:
- Black Red White Blue Yellow
 - Pink ¹ Green ¹ Purple ¹ Orange ¹
- ¹Available only with Spinergy LX
- Handrim Tabs:
- Available only with Aluminum Anodized and Plastic Coated Handrims*
- Short Tab _____ \$25
 - Long Tab _____ N/C

- Rear Tire Type**
- Full Profile**
- Street Pneumatic 1 3/8" E2211, E2212 _____ Std
Not available on 25" wheel size
 - Full Poly 1 3/8" E2218 _____ \$55
Not available on Superlight wheel or 25" or 26" wheel sizes
 - Pneumatic with Airless Insert 1 3/8" E2211, E2213 _____ \$100
Not available on Superlight wheel or 25" wheel size
 - Schwalbe HS-127 1 3/8" E2211, E2212 _____ \$125
Available only on 24" wheel size. Not available on Mag wheel.
 - Euro Trek Knobby 1 3/8" E2211, E2212 _____ \$95
Available only on 24" wheel size
 - Mountain 2" Black E2211, E2212 _____ \$95
Available only on 22", 24" and 25" MAXX Performance and Spinergy wheels
- Low Profile**
- Low Poly 1" E2218 _____ N/C
Not available on Superlight wheel, 25" or 26" wheel sizes or 350 lb Weight Capacity.
 - Sentinel Solid Tire 1" E2218 _____ \$95
Not available on Superlight wheel or 20" and 22" Mag wheels
 - High Pressure 1" E2211, E2212 _____ \$60
Not available on Mag wheel
 - Schwalbe One 1" E2211, E2212 _____ \$125
Available only on 24" and 25" wheel sizes. Not available on Mag wheel.
 - Schwalbe Marathon Plus Evolution 1" E2211, E2212 _____ \$140
Not available on Mag wheel, 20" or 22" wheel size

Rear Seat Height and Wheels continued

Handrim

Standard Handrims

- Aluminum Anodized _____ Std
- Superlight Aluminum Anodized _____ \$50
Available only on 24" and 25" wheel sizes. Available only on MAXX Performance, Superlight or Spinergy wheels. Not available with 350 lb Weight Capacity.
- Plastic Coated _____ \$95
- Projection (8 Vertical Projections) _____ \$325
Available only with full profile tires in combination with 22"-26" MAXX Performance or 24" Mag wheel
- Spinergy Flex Rim _____ \$600
Available only on 24" and 25" Spinergy LX wheels
- Omit Handrim _____ N/C

Natural Fit Handrims

Not available on 20" wheel size or 22" Mag wheel. Not recommended with Mag wheel.

- Natural Fit LT Regular Grip _____ \$465
- Natural Fit Standard Regular Grip _____ \$445
- Natural Fit LT Super Grip _____ \$515
- Natural Fit Standard Super Grip _____ \$495
- Natural Fit LT Omit Grip _____ \$395
- Natural Fit Standard Omit Grip _____ \$375

Handrim Option

- 3M Micro Grip _____ \$65
Not available on Plastic Coated or Projection Handrims

Wheel Locks

Wheel Lock

- Push to Lock _____ Std
- Push to Lock Flush Mount _____ N/C
- Pull to Lock _____ N/C
- Under Seat Scissor Lock _____ N/C
- Low Profile Scissor Lock _____ \$45
Not available with Low Poly, Full Poly or Sentinel Solid tires
- Short Thro Scissor Lock _____ \$110
- Omit Wheel Locks _____ N/C

Armrests

Armrest

- Height Adjustable T-Arm E0973
1.5" minimum Rear Wheel Spacing
 - Desk (8.5" - 13.5") _____ \$325
 - Full (8.5" - 13.5") _____ \$325
- Height Adjustable T-Arm Low E0973
1.5" minimum Rear Wheel Spacing
 - Desk (7.5" - 11.5") _____ \$325
 - Full (7.5" - 11.5") _____ \$325

Axle

- Quick Release _____ Std
- Quad Quick Release _____ \$185
Not available with Mag or Spinergy wheels

Center of Gravity Preset

Measured from front of back post to center of rear axle. 0" indicates axle will be directly under back post. See reference J on page 6.

Amputee Positions

<input type="checkbox"/> -1.5"	<input type="checkbox"/> -1.25"	<input type="checkbox"/> -1"	<input type="checkbox"/> -.75"
<input type="checkbox"/> -.5"	<input type="checkbox"/> -.25"	<input type="checkbox"/> 0"	<input type="checkbox"/> .25"

Standard Positions

<input type="checkbox"/> .5"	<input type="checkbox"/> .75"	<input type="checkbox"/> 1"	<input type="checkbox"/> 1.25"
<input type="checkbox"/> 1.5"	<input type="checkbox"/> 1.75"	<input type="checkbox"/> 2"	<input type="checkbox"/> 2.25"
<input type="checkbox"/> 2.5"	<input type="checkbox"/> 2.75"	<input type="checkbox"/> 3"	<input type="checkbox"/> 3.25"
<input type="checkbox"/> 3.5"	<input type="checkbox"/> 3.75"	<input type="checkbox"/> 4"	<input type="checkbox"/> 4.25"

Rear Wheel Spacing

Measured from outside of back to inside of rear tire. Adjustable out 1/2" from setting. 1.5" minimum Rear Wheel Spacing with Height Adjustable T-Arm Armrest. 1" minimum Rear Wheel Spacing with 0° Rear Wheel Camber. See reference K on page 6.

- .75" 1" 1.25" 1.5" 1.75"

Rear Wheel Camber

- 0° 2° 4° 6° 8°

Wheel Lock Options

- Extension Handles E0961 _____ \$77
Available only with Push to Lock or Pull to Lock Wheel Locks
- Grade Aids E0974 _____ \$150
Available only with Push to Lock Wheel Locks and Pneumatic or Pneumatic with Airless Insert tires

Side Guards

Not available with T-Arms or Swing Away with Side Guard Armrests

Direct Mount	
<input type="checkbox"/> Composite Side Guards K0108	\$150
Aluminum Straight K0108	
<input type="checkbox"/> Black Matte	\$195
<input type="checkbox"/> Match Frame Color	\$245
<i>Not available with Hydrographics</i>	
Aluminum Adjustable Fender with Foam Pad K0108	
<i>Not available with Adjustable Depth greater than 1/2"</i>	
<input type="checkbox"/> Black Matte	\$275
<input type="checkbox"/> Match Frame Color	\$325
<i>Not available with Hydrographics</i>	
<input type="checkbox"/> Carbon Fiber Adjustable Fender K0108	\$485
<i>Not available with Adjustable Depth greater than 1/2"</i>	

Accessories

Anti-Tip	
<input type="checkbox"/> Rear Anti-Tips E0971	\$145
<input type="checkbox"/> Omit Rear Anti-Tips	N/C
Calf Strap	
<input type="checkbox"/> Velcro Adjustable K0038	\$42
<input type="checkbox"/> Padded Velcro Adjustable K0038	\$68
<input type="checkbox"/> Bodypoint Padded AeroMesh K0038	\$118
Pelvic Positioning Belt	
<input type="checkbox"/> 1 1/2" Auto Buckle E0978	\$80
<input type="checkbox"/> 1 1/2" Padded Auto Buckle E0978	\$170
<input type="checkbox"/> 1" Auto Buckle E0978	\$80
<input type="checkbox"/> 1" Padded Auto Buckle E0978	\$170
<input type="checkbox"/> 2" Airline Buckle E0978	\$110
Bodypoint Evoflex - Pelvic Stabilizer E0978	\$188
Size: <input type="checkbox"/> Small (16-24") <input type="checkbox"/> Medium (19-29") <input type="checkbox"/> Large (25-37")	

Frame Color

Gloss Finish	
<input type="checkbox"/> Black	N/C
<input type="checkbox"/> Blue Skies	N/C
<input type="checkbox"/> Bubblegum Pink	N/C
<input type="checkbox"/> Candy Blue	N/C
<input type="checkbox"/> Candy Purple	N/C
<input type="checkbox"/> Candy Red	N/C
<input type="checkbox"/> Copper	N/C
<input type="checkbox"/> Granny Smith	N/C
<input type="checkbox"/> Metallic Orange	N/C
<input type="checkbox"/> Silver	N/C
<input type="checkbox"/> Yellow	N/C

Clamp-On Mount	
<input type="checkbox"/> Composite Side Guards K0108	\$175
Aluminum Straight K0108	
<input type="checkbox"/> Black Matte	\$225
<input type="checkbox"/> Match Frame Color	\$275
<i>Not available with Hydrographics</i>	
Aluminum Adjustable Fender with Foam Pad K0108	
<input type="checkbox"/> Black Matte	\$305
<input type="checkbox"/> Match Frame Color	\$355
<i>Not available with Hydrographics</i>	
<input type="checkbox"/> Carbon Fiber Adjustable Fender K0108	\$515

Fabric Side Guard	
<i>Not available when omitting Back Upholstery</i>	
<input type="checkbox"/> Fabric Side Guards K0108	\$85

Belt Mounting Kit	
<input type="checkbox"/> Ki Mobility Belt Mounting Clamp	\$35
<input type="checkbox"/> Bodypoint Belt Mounting Bracket Kit	\$108

Other Accessories	
<input type="checkbox"/> Tool Kit	\$45
<input type="checkbox"/> Removable Underseat Pouch	\$50
<input type="checkbox"/> Ki Mobility Backpack	\$65
<input type="checkbox"/> Luggage Carrier	\$85
<input type="checkbox"/> Black Neoprene Impact Guards	\$90
<input type="checkbox"/> Under Chair Carrier and Removable Bag	\$85
<i>Not available with scissor locks</i>	
<input type="checkbox"/> Clear Spoke Guards K0065	\$115
<i>Not available with Mag wheel</i>	
<input type="checkbox"/> Black Spoke Guards K0065	\$115
<i>Not available with 20" wheel size or Mag wheel</i>	
<input type="checkbox"/> Cane and Crutch Holder	\$195

Matte Finish	
<input type="checkbox"/> Charcoal Gray	N/C
<input type="checkbox"/> Emerald Green	N/C
<input type="checkbox"/> Iridescent Burgundy	N/C
<input type="checkbox"/> Iridescent Navy	N/C
<input type="checkbox"/> Military Green	N/C
<input type="checkbox"/> Retro Red	N/C
<input type="checkbox"/> Royal Blue	N/C
<input type="checkbox"/> Shadow Black	N/C
<input type="checkbox"/> Smooth White	N/C

Hydrographics	
Hydrographics	\$595
<input type="checkbox"/> AOR1 (Desert Digital Camo)	<input type="checkbox"/> Pink Camo
<input type="checkbox"/> Muddy Girl	<input type="checkbox"/> NEXT Camo
<input type="checkbox"/> Carbon Fiber ¹	<input type="checkbox"/> Zebra ¹

¹Select Hydrographics Base Color. Zebra not available in Black.

Hydrographics Base Color	N/C
<input type="checkbox"/> Black <input type="checkbox"/> Blue <input type="checkbox"/> Granny Smith <input type="checkbox"/> Pink <input type="checkbox"/> Purple <input type="checkbox"/> Red	

References

A. Seat Width See Diagram 1

Measured from outside of frame tube on one side to outside of frame tube on other side.

B. Seat Depth See Diagram 3

Measured from front of back posts to front edge of seat sling. Seat sling starts at beginning of bend at front of frame.

C. Sling Position See Diagram 2

In Performance Position (0°), sling is at front of frame to keep chair shorter and more maneuverable. The 1" or 2" seat sling position extends the frame in front of the sling by that amount. Adding more frame can improve stability and provide support to aid in transfers.

D. Front Frame Bend See Diagram 3

Angle between front seat tube and ground.

E. Footrest Taper See Diagram 1

Indicates bend of front frame creating footrest. Taper is measured from outside of seat frame to outside of front tube.

F. Front Seat Height See Diagram 1

Measured from floor to top of seat tube at front of seat sling.

G. Footrest Width See Diagram 1

Measured from inside of front tube to inside of front tube on other side. Footrest width is listed below for each seat width and taper.

Seat Width	Straight	1" Taper	2" Taper
Outside Measurement	2" Narrower than Seat Width	4" Narrower than Seat Width	6" Narrower than Seat Width
14"	11.5"	9.5"	
15"	12.5"	10.5"	
16"	13.5"	11.5"	9.5"
17"	14.5"	12.5"	10.5"
18"	15.5"	13.5"	11.5"
19"	16.5"	14.5"	12.5"
20"	17.5"	15.5"	13.5"
21"	18.5"	16.5"	14.5"
22"	19.5"	17.5"	15.5"

H. Seat to Footrest Length See Diagram 3

Measured from front edge of seat sling to top rear of footrest. Footrest length of at least 2.5" shorter than front seat height recommended.

I. Rear Seat Height See Diagram 3

Measured from floor to top of seat tube at front of back post. Rear seat height is custom to needs.

J. Center of Gravity Preset See Diagram 3

Measured from front of back post to center of rear axle. 0" indicates axle will be directly under back post.

K. Rear Wheel Spacing See Diagram 1

Measured from outside of seat back to inside of rear tire. Adjustable out .5" from setting.

Diagram 1




Diagram 2




Diagram 3

