



Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com

Company: _____

Ship To: _____

Account Number: _____ Date: _____

Address: _____

Quote PO Number: _____

City: _____ State: _____

RTS: _____

Zip: _____ Phone: _____

Marked For: _____

Email/Fax: _____

Select if chair quote has: Additional Ki Products with Different Quote/Order Number **Quote/Order Number:** _____

Effective 01/17/2022

Tsunami ALX K0005, E1235 **\$2070**
 250 lb weight capacity

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions. Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Frame

Seat Width

See reference A on page 6

- 14" _____ N/C
- 15" _____ N/C
- 16" _____ N/C
- 17" _____ N/C
- 18" _____ N/C
- 19" _____ \$125
- 20" E2201 _____ \$500
- 21" E2201 _____ \$500
- 22" E2201 _____ \$500

Seat Depth

See reference B on page 6

- 14" _____ N/C
- 15" _____ N/C
- 16" _____ N/C
- 17" _____ N/C
- 18" _____ N/C
- 19" _____ \$75
- 20" E2203 _____ \$550

Weight Capacity

- 250 lb _____ Std
- 350 lb (Heavy Duty) _____ \$750
Requires 6" casters

Transit

- Transit System _____ \$320
125 lb maximum weight capacity
- Omit Transit _____ N/C

Backrest

Back Post / Height Range / Back Height

Back Post	Price	Height Range	Back Height				
Height Adjustable	N/C	<input type="checkbox"/> Short	<input type="checkbox"/> 11 ²	<input type="checkbox"/> 12	<input type="checkbox"/> 13	<input type="checkbox"/> 14	
		<input type="checkbox"/> Medium	<input type="checkbox"/> 13	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17
		<input type="checkbox"/> Tall	<input type="checkbox"/> 17	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20	
Fixed Height ¹	\$85	<input type="checkbox"/> Fixed	<input type="checkbox"/> 9 ²	<input type="checkbox"/> 10 ²	<input type="checkbox"/> 11 ²	<input type="checkbox"/> 12 ³	<input type="checkbox"/> 13
			<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18

¹ Not available with back upholstery
² Not available with Fold Down Push Handle
³ Fold Down Push Handle requires Height Adjustable Rigidizer Bar

Sling Position

In Performance Position (0"), sling is at front of frame. The 1" or 2" sling position extends the frame in front of the sling by that amount. See reference C on page 6.

- 0" (Performance) _____ N/C
- 1" _____ N/C
- 2" (Classic) _____ N/C

Adjustable Depth

Adds the capability to grow seat in depth

- No Rear Adjustment _____ N/C
- 1/2" _____ N/C
- 1" _____ \$85
- 2" _____ \$165
- 3" _____ \$165

Long Sling

Available only with 1", 2" and 3" Adjustable Depth

- Long Sling _____ N/C
Sling will be maximum 2" longer than seat depth

Seat Option

- Standard Seat Sling _____ Std
- Solid Seat Pan E2231 _____ \$265
Seat Pan is 1" longer when Long Sling option is selected

Front Frame Bend

Angle between front seat tube and ground. See reference D on page 6.

- 90° _____ \$125
- 85° _____ N/C
- 80° _____ N/C
- 75° _____ N/C
- 70° _____ \$125

Push Handle

- Omit Push Handle _____ N/C
- Standard Push Handle _____ \$85
- Low Profile Push Handle _____ \$150
2" shorter than standard push handles
- Bolt-On Push Handle _____ \$125
Handles will be 1" - 2" above back height
- Fold Down Push Handle _____ \$100
Not available with 350 lb Weight Capacity

DCN0023.32 Pricing and specifications subject to change without notice

Backrest continued

- Rigidizer Bar**
- Adjustable Height Rigidizer Bar _____ \$75
 - Non-Adjustable Height Rigidizer Bar _____ N/C

- Locking Backrest**
- Locking Fold Down Back _____ \$95
Rogue style locking back. Not available with 350 lb Weight Capacity

- Seat Back Angle**
- Standard is 90° (perpendicular to floor). Minus (-°) closes back angle. Plus (+°) opens back angle.*
- Standard +3° +6° +9°
 - 3° -6° -9° -12°

Back Upholstery

- Available only with Height Adjustable Back Posts. If selecting both upholstery and Axiom Back, the upholstery will be installed.*
- Padded Velcro Adjustable
Not available with 350 lb Weight Capacity
- Standard Cover _____ Std
 - Air Mesh Breathable Cover _____ \$75
- Tension Adjustable E2611
- Standard Cover _____ \$320
 - Air Mesh Breathable Cover _____ \$390
- Omit Back Upholstery _____ - \$10

Military Patch

- US Military Service Patch _____ \$90
- Army Navy Marines Air Force Coast Guard
- US Purple Heart Patch _____ \$90
- Purple Heart with Military Service Patch Combo _____ \$145
- Army Navy Marines Air Force Coast Guard

Front Seat Height and Casters

Caster Size / Footrest Taper / Front Seat Height

Select Front Seat Height and Footrest Taper in row corresponding to Caster Size. 6" casters not available with 85° and 90° Front Frame Bend. Footrest Taper indicates bend of front frame creating footrest width. Taper is measured from outside of seat frame to outside of front tube. 2" Taper not available with 75° Front Frame Bend. Taper is not available with 70° Front Frame Bend. See reference E on page 6. Front Seat Height is measured from floor to top of seat tube at front of seat sling. Adjustable up or down by 1/2". See reference F on page 6.

Footrest Taper	Caster Size	Front Seat Height														
<input type="checkbox"/> Straight	3"	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5
	4"	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	
	5"	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5		
	6"	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5			
<input type="checkbox"/> 1" Taper	3"	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5						
	4"	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5							
	5"	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5								
	6"	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5									
<input type="checkbox"/> 2" Taper	3"	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5						
	4"	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5							
	5"	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5								
	6"	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5									

Caster

Caster Type	Caster Size	Caster Size			
		Price	3"	4"	5"
.75" Rollerblade	N/C	<input type="checkbox"/>	<input type="checkbox"/>		
.75" Lighted Rollerblade	\$70	<input type="checkbox"/>	<input type="checkbox"/>		
1" Poly	N/C		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
1" Poly Aluminum	\$110		<input type="checkbox"/>	<input type="checkbox"/>	
1" Pneumatic	\$80				<input type="checkbox"/>
1.5" Poly	\$55		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
1.5" Soft Roll Aluminum	\$185	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Fork Options

- Standard Fork _____ Std
- Slotted Fork _____ \$75
- Single Sided Fork _____ \$170
Not available with 6" Caster Size
- Frog Legs Suspension Fork E1015 _____ \$450
Not available with 350 lb Weight Capacity
 - User Weight 80 - 120 lbs
 - User Weight 120 - 160 lbs
 - User Weight 160 - 220 lbs
 - User Weight 220 - 250 lbs
- Frog Legs II Suspension System _____ \$450
Lightweight shock absorbing fork. Not available with 350 lb Weight Capacity.
- Caster Pin Locks K0073 _____ \$185
Available only with Slotted Fork

Footrest

Footplate

- Tubular "Open" Footplate _____ Std
Not available w/footrest widths greater than 17.5". See Reference G on page 6.
- Tubular Footplate with ABS Cover _____ \$30
Not available w/footrest widths greater than 15.5". See Reference G on page 6.
- Aluminum Angle Adjustable K0040 _____ \$95
- Aluminum Angle Adjustable with Riser K0040, K0037 _____ \$145
Not available with 350 lb Weight Capacity
- High Mount Aluminum Angle Adjustable K0040, K0037 _____ \$95
Not available with 350 lb Weight Capacity
- Angle Adjustable Flip Under K0040, K0037 _____ \$185
Not available with footrest widths less than 9.5". Not available with 350 lb Weight Capacity
- High Mount Angle Adjustable Flip Under K0040, K0037 _____ \$270
Not available with footrest widths greater than 15.5" or 350 lb Weight Capacity. See Reference G on page 6.

Footrest Protector

- Footrest Protector _____ \$25
Not available with Angle Adjustable Flip Under, High Mount Angle Adjustable Flip Under or High Mount Aluminum Angle Adjustable Footplates

Seat to Footrest Length

Measured from front edge of seat sling to top, rear of footrest. Footrest length at least 2.5" shorter than Front Seat Height recommended. See reference H on page 6.

<input type="checkbox"/> 8"	<input type="checkbox"/> 8.5"	<input type="checkbox"/> 9"	<input type="checkbox"/> 9.5"	<input type="checkbox"/> 10"	<input type="checkbox"/> 10.5"
<input type="checkbox"/> 11"	<input type="checkbox"/> 11.5"	<input type="checkbox"/> 12"	<input type="checkbox"/> 12.5"	<input type="checkbox"/> 13"	<input type="checkbox"/> 13.5"
<input type="checkbox"/> 14"	<input type="checkbox"/> 14.5"	<input type="checkbox"/> 15"	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"
<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"

Rear Seat Height and Wheels

Wheel Size / Rear Seat Height

Select Rear Seat Height in row corresponding to desired Wheel Size. For proper setup, selection is required when omitting rear wheels. Measured from floor to top of seat tube at front of back post. Rear Seat Height is custom to needs. See reference I on page 6.

Wheel Size	Rear Seat Height													
22"	<input type="checkbox"/> 13.5 ¹	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20
24"	<input type="checkbox"/> 14.5 ¹	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20		
25"	<input type="checkbox"/> 15 ¹	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20			
26"	<input type="checkbox"/> 15.5 ¹	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20				

¹ Available only with Low Profile tires

Camber Bar Tower

- Standard Tower _____ N/C
- Mx Tower _____ \$100

Standard Tower allows rear seat height adjustment up or down in 1/4" increments. 1" adjustment available; 1/2" in lowest end of range. Mx Tower required with select low seat heights or when large rear seat height adjustment range is required. Provides 3" adjustment range in 1/4" increments.

Rear Wheel

		Rear Wheel Size				
		Price	22"	24"	25"	26"
Rear Wheel Type	MAXX Performance Spoke	\$125	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	MAXX Performance Spoke	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Mag	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Spinergy Spox	\$845	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Spinergy Options

Select spoke color and handrim tab for Spinergy wheels

Spoke Color:

- Black Red White Blue Yellow

Handrim Tab:

Available only with Aluminum Anodized and Plastic Coated Handrims

- Short Tab _____ \$25
- Long Tab _____ N/C

Rear Tire Type

Full Profile

- Street Pneumatic 1 3/8" E2211, E2212 _____ Std
Not available on 25" wheel size
- Full Poly 1 3/8" E2218 _____ \$55
Not available on 25" or 26" wheel sizes
- Pneumatic with Airless Insert 1 3/8" E2211, E2213 _____ \$85
Not available on 25" wheel size
- Schwalbe HS-127 1 3/8" E2211, E2212 _____ \$125
Available only on 24" wheel size. Not available on Mag wheel.
- Euro Trek Knobby 1 3/8" E2211, E2212 _____ \$95
Available only on 24" wheel size
- Mountain 2" Black E2211, E2212 _____ \$125
Available only on 22", 24" and 25" MAXX Performance and Spinergy wheels

Low Profile

- Low Poly 1" E2218 _____ N/C
Not available on 25" or 26" wheel sizes or 350 lb Weight Capacity
- Sentinel Solid Tire 1" E2218 _____ \$95
Not available on 22" Mag wheel
- High Pressure 1" E2211, E2212 _____ \$80
Not available on Mag wheel
- Schwalbe Marathon Plus Evolution 1" E2211, E2212 _____ \$140
Not available on Mag wheel or 22" wheel size

Rear Seat Height and Wheels continued

Handrim

Standard Handrims

- Aluminum Anodized _____ Std
- Plastic Coated _____ \$112
- Projection (8 Vertical Projections) _____ \$325
Available only with full profile tires in combination with MAXX Performance or 24" Mag wheel
- Omit Handrim _____ N/C

Natural Fit Handrims

Not available on 20" wheel size or 22" Mag wheel. Not recommended with Mag wheel.

- Natural Fit LT Regular Grip _____ \$465
- Natural Fit Standard Regular Grip _____ \$445
- Natural Fit LT Super Grip _____ \$515
- Natural Fit Standard Super Grip _____ \$495
- Natural Fit LT Omit Grip _____ \$395
- Natural Fit Standard Omit Grip _____ \$375

Handrim Option

- 3M Micro Grip _____ \$65
Not available on Plastic Coated or Projection Handrims

Axle

- Quick Release _____ Std
- Quad Quick Release _____ \$185
Not available with Mag or Spinergy wheels

Wheel Locks

Wheel Lock

- Push to Lock _____ Std
- Push to Lock Flush Mount _____ N/C
- Pull to Lock _____ N/C
- Under Seat Scissor Lock _____ N/C
- Low Profile Scissor Lock _____ \$65
Not available with Low Poly, Full Poly or Sentinel Solid tires
- Short Thro Scissor Lock _____ \$110
- Omit Wheel Locks _____ N/C

Armrests

Armrest

- Height Adjustable T-Arm (8.5" - 13.5") E0973
1.5" minimum Rear Wheel Spacing
- Desk _____ \$325
 - Full _____ \$325
- Height Adjustable T-Arm Low (7.5" - 11.5") E0973
1.5" minimum Rear Wheel Spacing
- Desk _____ \$325
 - Full _____ \$325

Center of Gravity Preset

Measured from front of back post to center of rear axle. 0" indicates axle will be directly under back post. See reference J on page 6.

Amputee Positions

<input type="checkbox"/> -1.5"	<input type="checkbox"/> -1.25"	<input type="checkbox"/> -1"	<input type="checkbox"/> -.75"
<input type="checkbox"/> -.5"	<input type="checkbox"/> -.25"	<input type="checkbox"/> 0"	<input type="checkbox"/> .25"

Standard Positions

<input type="checkbox"/> .5"	<input type="checkbox"/> .75"	<input type="checkbox"/> 1"	<input type="checkbox"/> 1.25"
<input type="checkbox"/> 1.5"	<input type="checkbox"/> 1.75"	<input type="checkbox"/> 2"	<input type="checkbox"/> 2.25"
<input type="checkbox"/> 2.5"	<input type="checkbox"/> 2.75"	<input type="checkbox"/> 3"	<input type="checkbox"/> 3.25"
<input type="checkbox"/> 3.5"	<input type="checkbox"/> 3.75"	<input type="checkbox"/> 4"	<input type="checkbox"/> 4.25"

Rear Wheel Spacing

Measured from outside of back to inside of rear tire. Adjustable out 1/2" from setting. 1.5" minimum Rear Wheel Spacing with Height Adjustable T-Arm Armrest. 1" minimum Rear Wheel Spacing with 0° Rear Wheel Camber. See reference K on page 6.

- .75" 1" 1.25" 1.5" 1.75"

Rear Wheel Camber

- 0° 2° 4° 6° 8°

Wheel Lock Options

- Extension Handles E0961 _____ \$77
Available only with Push to Lock or Pull to Lock Wheel Locks
- Grade Aids E0974 _____ \$150
Available only with Push to Lock Wheel Locks and Pneumatic or Pneumatic with Airless Insert tires

Side Guards

Not available with T-Arms or Swing Away with Side Guard Armrests

Direct Mount

- Composite Side Guards - Direct Mount Receiver K0108 ___ \$150
- Aluminum Straight K0108
 - Black Matte _____ \$195
 - Match Frame Color _____ \$245
- Aluminum Adjustable Fender with Foam Pad K0108
Not available with Adjustable Depth greater than 1/2"
 - Black Matte _____ \$275
 - Match Frame Color _____ \$325
- Carbon Fiber Adjustable Fender K0108 _____ \$485
Not available with Adjustable Depth greater than 1/2"

Clamp-On Mount

- Composite Side Guards - Clamp-on Receiver K0108 _____ \$175
- Aluminum Straight K0108
 - Black Matte _____ \$225
 - Match Frame Color _____ \$275
- Aluminum Adjustable Fender with Foam Pad K0108
 - Black Matte _____ \$305
 - Match Frame Color _____ \$355
- Carbon Fiber Adjustable Fender K0108 _____ \$515

Fabric Side Guard

Not available when omitting Back Upholstery

- Fabric Side Guards K0108 _____ \$95

Accessories

Anti-Tip

- Rear Anti-Tips E0971 _____ \$150
- Omit Rear Anti-Tips _____ N/C

Calf Strap

- Velcro Adjustable K0038 _____ \$42
- Padded Velcro Adjustable K0038 _____ \$68
- Bodypoint Padded AeroMesh K0038 _____ \$118

Pelvic Positioning Belt

- 1 1/2" Auto Buckle E0978 _____ \$80
- 1 1/2" Padded Auto Buckle E0978 _____ \$170
- 2" Airline Buckle E0978 _____ \$110
- Bodypoint Evoflex - Pelvic Stabilizer E0978 _____ \$188
 Size: Small (16-24") Medium (19-29") Large (25-37")

Belt Mounting Kit

- Ki Mobility Belt Mounting Clamp _____ \$35
- Bodypoint Belt Mounting Bracket Kit _____ \$108

Other Accessories

- Tool Kit _____ \$45
- Removable Underseat Pouch _____ \$50
- Ki Mobility Backpack _____ \$65
- Luggage Carrier _____ \$85
- Black Neoprene Impact Guards _____ \$90
- Under Chair Carrier and Removable Bag _____ \$85
Not available with scissor locks
- Clear Spoke Guards K0065 _____ \$115
Not available with Mag wheel
- Black Spoke Guards K0065 _____ \$115
Not available with 20" wheel size or Mag wheel
- Cane and Crutch Holder _____ \$195

Frame Color

Gloss Finish

- Black _____ N/C
- Blue Skies _____ N/C
- Bubblegum Pink _____ N/C
- Candy Blue _____ N/C
- Candy Purple _____ N/C
- Candy Red _____ N/C
- Copper _____ N/C
- Granny Smith _____ N/C
- Metallic Orange _____ N/C
- Silver _____ N/C
- Yellow _____ N/C

Matte Finish

- Charcoal Gray _____ N/C
- Emerald Green _____ N/C
- Iridescent Burgundy _____ N/C
- Iridescent Navy _____ N/C
- Military Green _____ N/C
- Retro Red _____ N/C
- Royal Blue _____ N/C
- Shadow Black _____ N/C
- Smooth White _____ N/C

References

A. Seat Width See Diagram 1

Measured from outside of frame tube on one side to outside of frame tube on other side.

B. Seat Depth See Diagram 3

Measured from front of back posts to front edge of seat sling. Seat sling starts at beginning of bend at front of frame.

C. Sling Position See Diagram 2

In Performance Position (0°), sling is at front of frame to keep chair shorter and more maneuverable. The 1" or 2" seat sling position extends the frame in front of the sling by that amount. Adding more frame can improve stability and provide support to aid in transfers.

D. Front Frame Bend See Diagram 3

Angle between front seat tube and ground.

E. Footrest Taper See Diagram 1

Indicates bend of front frame creating footrest. Taper is measured from outside of seat frame to outside of front tube.

F. Front Seat Height See Diagram 1

Measured from floor to top of seat tube at front of seat sling.

G. Footrest Width See Diagram 1

Measured from inside of front tube to inside of front tube on other side. Footrest width is listed below for each seat width and taper.

Seat Width	Straight	1" Taper	2" Taper
Outside Measurement	2" Narrower than Seat Width	4" Narrower than Seat Width	6" Narrower than Seat Width
14"	11.5"	9.5"	
15"	12.5"	10.5"	
16"	13.5"	11.5"	9.5"
17"	14.5"	12.5"	10.5"
18"	15.5"	13.5"	11.5"
19"	16.5"	14.5"	12.5"
20"	17.5"	15.5"	13.5"
21"	18.5"	16.5"	14.5"
22"	19.5"	17.5"	15.5"

H. Seat to Footrest Length See Diagram 3

Measured from front edge of seat sling to top rear of footrest. Footrest length of at least 2.5" shorter than front seat height recommended.

I. Rear Seat Height See Diagram 3

Measured from floor to top of seat tube at front of back post. Rear seat height is custom to needs.

J. Center of Gravity Preset See Diagram 3

Measured from front of back post to center of rear axle. 0" indicates axle will be directly under back post.

K. Rear Wheel Spacing See Diagram 1

Measured from outside of seat back to inside of rear tire. Adjustable out .5" from setting.

Diagram 1




Diagram 2




Diagram 3

